

Gynaecology Plinths




GY02E/GY04E

Engineered for modern clinical settings, the GP Gynaecological Plinth offers enhanced clinician access and patient comfort through its ergonomic design and electric height-and-tilt adjustment, including Trendelenburg positioning. The fully electric model features versatile leg support options, foot stirrups, knee troughs, or upholstered footpads and includes a sturdy stainless-steel tray for convenience and hygiene.

It also converts easily into a standard two-section plinth thanks to a removable leg section. Its low entry height, gas-spring assisted backrest, retractable castors, and adjustable foot for uneven floors make it highly adaptable in fast-paced care environments.

The Ultra Gynaecological Plinth represents the pinnacle of comfort and clinical versatility. Featuring three independent electric motors for controlling height, tilt, and backrest, it facilitates smooth Trendelenburg positioning. Standard fittings include adjustable knee troughs or foot stirrups, robust stainless-steel tray, and soft, durable double-thickness upholstery with infection- and fire-resistant properties. UK-manufactured and CE-marked.

Code	Mechanism	SWL	Surface	Footrest Adj	Backrest/Mid Section Adj
GY02ET	Electric	260kg	173 x 63 cm	N/a	85°
GY02ES	Electric	260kg	173 x 63 cm	N/a	85°
GY02EF	Electric	260kg	173 x 63 cm	N/a	85°
GY04ET	Electric	260kg	184 x 63 cm	N/a	85°
GY04ES	Electric	260kg	184 x 63 cm	N/a	85°
GY04EF	Electric	260kg	184 x 63 cm	N/a	85°




Accessories	Code
Troughs, Stirrups or Foot Pads	GY03/GY01/GY05
Breathing Hole	BR02
Paper Roll Holder	PR01
Adjustable Arms	DE01
Detachable Head Rest	DHR01
Electric Backrest	EB01
Fold Down Arm Loops (GY02 only)	AL01
Footswitch	FO02 or FO03 (model permitting)
Positioning Wedge Set	WEDGESET



14 Standard Colours



BEIGE



CITRUS GREEN




APPLE MINT



CANDY




RASPBERRY CRUSH



DAMSON



PILLAR BOX RED




WHITE



ASH GREY



GUN METAL GREY




BLACK



DARK BLUE



MID BLUE



SKY BLUE