

# Bariatric Plinths






## **BARIO275/100, BARIO375/100, BARIO375DE/100DE**

The Medi-Plinth Bariatric Plinth range is designed to meet the growing demand for high-capacity treatment equipment, offering robust, feature-rich solutions that ensure both patient comfort and clinical efficiency. Each model supports a maximum safe working load of 320 kg and incorporates a secure electric lift system powered by Linak “High Reliability” actuators. Optional extra-wide surfaces of 75 cm or 100 cm further enhance patient accessibility and comfort. Adjustable feet and an extra-low entry height facilitate seamless wheelchair transfers, and retractable castors allow easy movement between treatment areas.

The 2-Section Bariatric Couch delivers a compact yet supportive platform with electric backrest as standard. The 3-Section Bariatric Couch builds on this functionality with an additional section that enables greater flexibility in patient positioning, including elevated hip and leg adjustments for more complex procedures. The Drop-End Bariatric Couch provides the added benefit of transforming the plinth into a chair to support a wider variety of patient care routines and improve clinical accessibility. UK-manufactured and CE-marked.

Code	Mechanism	SWL	Surface	Footrest Adj	Backrest Adj
BARIO275	Electric	320kg	186 x 75 cm	N/a	-20° to +85°
BARIO2100	Electric	320kg	186 x 100 cm	N/a	-20° to +85°
BARIO375	Electric	320kg	186 x 75 cm	-20° to +85°	-20° to +85°
BARIO3100	Electric	320kg	186 x 100 cm	-20° to +85°	-20° to +85°
BARIO3100DE	Electric	320kg	186 x 100 cm	-85°	-15° to +70°
BARIO375DE	Electric	320kg	186 x 100 cm	-85°	-15° to +70°



Accessories	Code
Footswitch	FO02
Extra Wide Upholstery- large format couches only (120cm)	EW3
Breathing Hole	BR02
Detachable Head Rest (75cm wide)	DHR02
Detachable Head Rest (100cm wide)	DHR03
Bariatric Armrests	BE01
Electric Foot	+EB03



Find our products on the NHS Supply Chain Catalogue

## 14 Standard Colours



BEIGE



CITRUS GREEN



APPLE MINT



CANDY




RASPBERRY CRUSH



DAMSON



PILLAR BOX RED




WHITE



ASH GREY



GUN METAL GREY




BLACK



DARK BLUE



MID BLUE



SKY BLUE