

Tilt Tables



T01E, T02E, T03E


Medi-Plinth's Tilt Tables are engineered for reliable patient positioning in clinical environments, featuring dual-motor electric systems for smooth height and tilt adjustments.

Standard Tilt Table- This versatile model features large lockable castors for patient transport, integrated patient handles, an inclinometer, and three Velcro harnesses. A fixed footboard supports comfort and stability during tilting.

2-Section Tilt Table- Building on the standard model, the table includes an adjustable backrest, allowing the unit to function both as a tilt table and a general examination couch. Additional features include a fixed footboard, Velcro harnesses and patient handles.

Split-Leg Tilt Table- This variant incorporates a split-leg section that enables separate strapping of each leg for enhanced patient security—especially beneficial for those with limited stability. It retains the full tilt and height control, large transport castors, passenger handles, inclinometer, and protective harnesses seen in the standard version.

Code	Mechanism	SWL	Surface	Footrest Adj	Backrest Adj
T02E	Electric	180kg	186 x 63 cm	N/a	N/a
T01E	Electric	180kg	186 x 63 cm	N/a	N/a
T03E	Electric	180kg	186 x 63 cm	N/a	-15° to 85°




Accessories	Code
Footswitch	FO02
Extra Wide Upholstery (75cm)	EW1
Breathing Hole	BR01
Electric Backrest (2 Section Only)	EB01
Tilting Foot Board	TF01
Tilting Split Foot Board	TF02
Tilting Patient Table	PT02
Velcro Harness (Set of 3 for Standard Tilt Tables)	VE03
Velcro Harness (Set of 3 for Split Leg Tilt Tables)	VE04



14 Standard Colours



BEIGE



CITRUS GREEN




APPLE MINT



CANDY




RASPBERRY CRUSH



DAMSON



PILLAR BOX RED




WHITE



ASH GREY



GUN METAL GREY




BLACK



DARK BLUE



MID BLUE



SKY BLUE