

Neonatal Chair



Neonatal Chair 87212

The Neonatal Chair is thoughtfully engineered to support both parents and caregivers during extended newborn care in neonatal units and maternity settings. It offers electric height adjustability from 47 to 85 cm, facilitating optimal positioning for breastfeeding and skin-to-skin contact.

Compact yet functional, the chair's ergonomic design takes up significantly less space than a hospital bed without compromising usability.

It provides exceptional back and neck support through anatomically contoured upholstery, and is built for clinical durability with edge-banded backrests to prevent dirt accumulation and simplify cleaning.

Safety and compliance are paramount: the Neonatal Chair is CE marked under MDR 2017/745 and meets the latest EMC (electromagnetic compatibility) standards for medical electrical devices



Code	Headrest	Backrest (W x H)	Seat (W x D)	Footrest
87212+COLOUR CODE	52 x 35cm	52 x 52cm	52 x 48cm	44 x 33cm

Code	Armrest Pad (W x L)	Seat Height	External Dims	SWL
87212+COLOUR CODE	12 x 63cm	44-84cm	180 x 76cm	140kg

Colour Code

