

# VELA Mammography Chair



PUSH BAR MAKES IT EASY  
TO MOVE THE CHAIR AND  
POSITION THE PATIENT  
CORRECTLY

ARMRESTS MAKE IT EASY  
FOR THE PATIENT TO SIT  
DOWN AND GET UP FROM  
THE CHAIR

ELECTRIC HEIGHT  
ADJUSTMENT ENSURES  
FAST AND CORRECT POSI-  
TIONING OF THE PATIENT  
RELATIVE TO THE MACHINE

THE CHAIR IS EASILY LOCKED WITH  
A FOOT BRAKE, WHICH ENSURES A  
STATIONARY AND STABLE PATIENT

FOOTREST PROVIDES THE  
PATIENT WITH INCREASED  
STABILITY REGARDLESS OF  
THE HEIGHT SETTING

BATTERY INDICATOR

Item no. 70032

The VELA Mammography Chair is optimal for patients who need to sit down during examinations, e.g. women prone to fainting, wheelchair users or patients with diseases affecting the musculoskeletal system.

- :: **Quick and secure patient positioning** for the examination device using the chair's smooth-running wheels, push bar and foot brake
- :: More **accurate examinations** for anxious or unsteady patients as they are more able to relax when seated
- :: **Optimal positioning** of the patient with electric height adjustment, which makes it **easy to vary the chair height** to suit the equipment and staff
- :: **Good working position for the staff** because regardless of the patient's height they can be seated at the right height for the staff, **reducing stress in the back, neck and shoulders**
- :: **Forward-leaning position** easily achieved and secured by adjusting the backrest
- :: The patient can be **transferred** between, for example a wheelchair and the patient chair due to **foldable armrests**