

Children's Rollators



Fixi & Fox Rollators

The Fixi and Fox paediatric rollators are purpose-built to provide safe, stable mobility support for children in clinical, rehabilitation, and home care settings.

Constructed from lightweight yet robust materials, both models ensure durability while remaining easy to manoeuvre. Their compact frames are designed for indoor and outdoor use, enabling young users to move with greater confidence and independence.

Height-adjustable handles allow precise fitting to each child's needs, while ergonomic grips enhance comfort and control during use.

The Fixi features a standard mobility frame for everyday walking support, whereas the Fox incorporates a forearm support frame for children requiring additional upper-body stability. Both models are equipped with reliable braking systems for safety, easy folding mechanisms for transport and storage, and smooth-rolling wheels that perform well across different floor types.



Code	Width	Depth	Height of Handle	Weight	Height Adj Steps	SWL
186.00.1*	68cm	58cm	61-70cm	6.2kg	5	100kg
186.00.2*	68cm	58cm	70-78cm	6.3kg	5	100kg

Configuration

186.00.10: Fixi Rollator for Children, Handle Height 61-70cm

186.00.15: Fixi Rollator for Children, Inc Seat, Handle Height 61-70cm

186.00.20: Fox Rollator for Children, Handle Height 70-78cm

186.00.25: Fox Rollator for Children, Inc Seat, Handle Height 70-78cm

