

# Flexline Anaesthesia Cart



## FLANES1XP Anaesthesia Cart

The Flexline Anaesthesia Cart is purpose-built for high-security clinical environments where controlled storage and seamless workflow are essential.

Crafted from lightweight yet robust polymer with integrated antimicrobial Microban protection, the cart ensures hygiene is maintained between cleanings while offering gentle maneuverability via built-in handles and smooth polyurethane castors

Security is a core feature with each self-closing, full-access drawer is equipped with an individual key lock and serialised plastic seals to uphold protocols. An optional electronic keypad with LCD display elevates access control, allowing enhanced security management.

Storage is versatile and efficient, combining a range of drawer depths (7.5cm, 15.2cm, and 22.9cm) with included dividers and accessories such as ampoule trays, side bins, a waste basket, peel-pouch catheter holders, overbridge hanger rails, and tilt bins. A spacious top surface supports procedural prep or documentation with ease.

Code	Width	Length	Height	Weight
FLANES1XP	55.6cm	81.9cm	114.9cm	91kg

## Configuration

FLANES1XP: Anaesthesia Cart with Key Lock System & Electronic Keypad Locking System, 1 x 75mm Key Locking Drawer, 2 x 7.5cm, 1 x 15.2cm & 1 x 22.9cm Deep Drawers with Blue Drawer Pulls. Accessories: 1 x 75mm Drawer Divider Kit, 2 x 152/229mm Drawer Divider Kit, 4 Locking Side Bins, Wastebasket & Holder, Peel Pouch Catheter Holder, Overbridge With 2 Hanger Rails & 1 Set of 4 Tilt Bins & 1 Set of 6 Tilt Bins Security Seals

Accessory	Code
3" Drawer Divider Kit for Flexline	FL113
6"/9" Drawer Divider Kit for Flexline and Lifeline Carts	FL116
Flexline Oxygen Tank Holder, All Size Carts, 4" Diameter Tank	FL234
Adjustable Defibrillator Tray for Lifeline and Flexline Carts	LEC304