

Geriatric Recliner Chair



6700 Hydraulic Recliner Chair

The Hydraulic Geriatric Recliner Chair delivers optimal comfort, functionality, and hygiene for clinical and care environments. Engineered with patient support in mind, the chair incorporates a high-density, easy-to-clean Future fabric upholstery that is waterproof, durable, resistant to disinfectants, and suitable even for outdoor use.

Constructed on an epoxy-coated steel frame that meets MDR and CE medical device certification standards. Manoeuvrability is ensured by Ø125mm locking castors, with the option to upgrade to larger Ø300mm front castors for challenging surfaces.

Ergonomic adjustments use reliable gas spring mechanisms to control backrest, seat tilt, leg section, and extendable footrest, complemented by fold-back armrests and neck cushions to enhance patient positioning.

Integrated practical features such as sturdy push handles and a detachable meal tray (adjustable in depth) make it suitable for multi-functional care roles. Accessories like IV pole adapters and seat belts further extend its adaptability within clinical workflows.

Code	Total Width	Seat Height	Seat Width	Backrest adj	Footrest length adj
6700	67cm	55cm	49cm	26°	14cm

Code	Leg Section adj	Headrest Height adj	SWL	Weight w/o Accessories
6700	-2 - 47°	9cm	150kg	59kg

Configuration

6700: Geriatric Recliner Chair, Hydraulic- inc Neck Cushion, Meal Tray, Fold Away Foot Rest, Arm Rests & Push Handles.

Accessories	Code
Ø300mm Wheels in Front	WLFW300/125GH
Seat Belt	6700150
IV Pole Adapter	GH20871
IV Pole, 4 Hooks. Height Adjustment with One Hand - Adapter Required	60121AMI