

Crutches



Safe-In Crutches Range

Rebotec's Safe-In forearm crutches deliver a blend of ergonomic design, durability, and safety, ideal for healthcare and rehabilitation settings.

The standard-size version features a near-closed forearm cuff with a hinge, reducing wrist fatigue and enabling a more natural arm swing. It includes an ergonomic handle with front and rear reflectors for improved visibility and ten-step height adjustability, supporting user weights up to 130 kg.

For paediatric use, the Safe-In Kids model offers the same supportive and adjustable design in a lightweight, colourful format tailored to young users' needs. It accommodates up to 100 kg and features three height adjustments, making it a practical solution for child rehab and mobility.

Rebotec also offers a Safe-In Extra Long version, which extends handle height from 85 to 107 cm to accommodate taller individuals, while retaining the same ergonomic and weight-bearing features.



Code	Height from Handle to Cuff	Height of Handle	Diameter Ø in mm	SWL
116.10.20	22-25cm	85-107cm	22/19	130kg
116.10	22-25cm	75-97cm	22/19	130kg
120.50.10	16-21cm	75-97cm	20/16.9	130kg
117.10	22-25cm	74-97cm	22/19	130kg

Configuration

116.10.20: Safe-In Crutch Size Extra Long

116.10: Safe-In Crutch Size Standard

120.50.10: Safe-In-Light Crutch Size Youth

117.10: Safe-In-Soft Crutch with Ergonomic Handle Size Standard



Find our products on the NHS Supply Chain Catalogue

