


SS Dressing Trolleys



DRT001/05/06/10/14/18/21

The Shuttleworth DRT range of stainless steel dressing trolleys has been engineered to meet the practical and hygiene demands of clinical environments, making them a reliable choice for NHS and private healthcare procurement.

Built from high-grade stainless steel, they offer a durable and corrosion-resistant solution that withstands daily use and rigorous cleaning protocols, ensuring long service life and compliance with infection control standards.

Each trolley is designed to support efficient workflow, with a stable welded frame and removable stainless steel trays that can be easily sterilised or replaced. Smooth-running swivel castors allow for effortless mobility, even in restricted spaces, while lockable options provide additional stability when in use.

Available in a range of sizes and configurations, the DRT trolleys deliver flexibility for general wards, treatment rooms, and theatres.

The range represents excellent value for procurement teams, combining robust construction with low maintenance requirements.

Code	Width	Depth	Height
DRT001	47cm	47cm	89.5cm
DRT005	46cm	46cm	89.6cm
DRT006	61cm	46cm	89cm
DRT010	47cm	47cm	89cm
DRT014	46cm	46cm	91.7cm
DRT018	122cm	46cm	91.5cm
DRT021	47cm	47cm	89.8cm

Configuration

DRT001: Stainless Steel Dressing Trolley, Removable Shelves, 3 Edges Up 1 Down

DRT005: Stainless Steel Dressing Trolley, Fixed Flat Shelves

DRT006: Stainless Steel Dressing Trolley, Fixed Flat Shelves

DRT010: Stainless Steel Dressing Trolley, Removable Shelves, Flat Shelves

DRT018: Stainless Steel Dressing Trolley, Fixed Shelves 3 Edges Up 1 Down


DRT021: Mild Steel Dressing Trolley, Removable Shelves, 3 Edges Up 1 Down

Accessories

Code

Four Rubber Buffers added to each corner

DRT020.1



Handling

The DRT Stainless Steel Dressing Trolleys are engineered for smooth and reliable performance in everyday clinical use.

Designed with a compact frame and easy-glide castors, the trolleys can be moved effortlessly between treatment rooms or wards, enabling efficient organisation of instruments and supplies.

The ergonomic push handle and lightweight stainless steel structure further enhance manoeuvrability while maintaining strength and durability in demanding healthcare settings.



Longevity

Built entirely from high-grade stainless steel, the DRT Trolleys provide exceptional durability and resistance to corrosion, ensuring a long service life even under frequent cleaning and disinfection.

Each model is designed with smooth, wipe-clean surfaces that support infection control compliance and simplify day-to-day maintenance. The reinforced shelving provides dependable support for medical equipment and consumables, while raised edges help to secure items during transport, reducing the risk of spills.



Safety

Safety and reliability are central to the DRT Dressing Trolleys. Heavy-duty castors with locking mechanisms ensure stability when stationary, giving clinicians confidence during delicate procedures or when accessing critical supplies.

The robust stainless steel frame supports consistent use in high-intensity care environments, and all materials comply with strict NHS infection control standards. The DRT Trolleys deliver consistent, safe performance for hospitals, clinics, and community healthcare settings.