

# manuthera<sup>®</sup> 242


Loved by professionals – Master of Spine Techniques



## Main features:

1. Dynamic patient surface for the most effective cervical, thoracic and lumbar treatment
2. Electronic height and tilt adjustment for both more efficient decompression treatment, ergonomics and easy access for patients
3. New improved patient comfort and safety
4. Stylish and modern new carbon-colored design

Manuthera 242 isn't just any table – it's a unique and versatile companion for most demanding professionals dedicated to manual therapy and patient care.

Technical data and features		Manuthera 242 dimensions
Thorax section elements	2	
Height adjustment	51–106 cm	
Weight	170 kg	
Tilt Up / Down	Electrical 52°/32°	
Tilt Up/Down function	Electronic	
Lower body rotation	Yes	
Lower body section movement (3D)	+21.5°/max -40° (vertical), 0..18° (horizontal)	
Upper thorax -section movement (3D)	+15°/-13° (vertical), 0..18° (horizontal)	
Total lock function for upper/lower body	Yes	
Castors	Yes	
Safe working load (SWL)	150 kg	
Lifting motors (push only)	2 (Motor push to lift, gravitational lowering)	
Drop functions	Traction drop	
Automatic anti-squeeze safety function	Yes	
Safety edge	Yes	
Adjustable legs and paws	Yes	
Safety lock functions	Automatic (2 h) or Manual / Manually operated	
Lying surface	Antimicrobial upholstery	
U-shaped face element and breathing hole cover	Yes	
Manipulation handle	Yes	
Body strap	Yes	
Country of Origin	Finland	
CE Marked	Yes The device is a CE Marked medical device according to MDR 2017/745 class 1)	
Manufacturer	Lojer Oy	

