

# Mobile Commode Chairs



## Bonn, Hamburg & Boston Commode Chairs

The Rebotec Bonn, Hamburg, and Boston mobile commode chairs are manufactured from high-strength, fully synthetic materials with seamless joints to prevent dirt accumulation and ensure optimal hygiene.

Designed for durability in high-use environments, the frames are completely rustproof and easy to disinfect, making them suitable for hospitals, care homes, and community care.

All models feature retractable armrests for easier lateral transfers, ergonomically contoured backrests for patient comfort, and integrated push handles for effortless manoeuvrability.

Four swivel castors with brakes provide secure positioning during use, while height-adjustable swing-away footrests further enhance patient support. The commodes are supplied with a removable commode pan and lid, ensuring safe and hygienic waste handling.

Fit with four 125mm swivel castors, each with an individual brake, delivering stability and smooth handling.

Code	Width	Depth	Seat Width	Seat Height	Seat Depth
343.**.00 (Bonn)	57cm	87cm	46cm	53cm	45cm
356.54.00 (Hamburg)	56cm	87cm	46cm	54-59cm	48cm
342.54.00 (Boston)	62cm	89cm	52cm	53cm	48cm

Code	Toilet Height w/out Bucket	Toilet Height w/ Bucket	Castors	SWL
343.**.00 (Bonn)	47cm	44cm	12.5cm (5")	130kg
356.54.00 (Hamburg)	45-50cm	48-53cm	12.5cm (5")	130kg
342.54.00 (Boston)	44cm	47cm	12.5cm (5")	150kg

## Configuration

343.54.00: Bonn Mobile Commode Chair, Fixed Height

343.54.22: Bonn Mobile Commode Chair with Hygiene Opening, Fixed Height

356.54.00: Hamburg Mobile Commode Chair- Adjustable Height

342.54.00: Boston Wide Mobile Commode Chair- Fixed Height