

Bo-Plinths




BO01E

The Medi-Plinth Two-Section Bobath Plinth is designed to meet the rigorous needs of neurologically complex and head-injury patients, combining strength, flexibility, and ergonomic treatment support. Its robust, wide patient surface provides enhanced clinician access and confidence during therapy and examination procedures.

Dual gas-spring assisted back section facilitates easy adjustment, while a heavy-duty adjustable foot helps maintain stability on uneven floors. The extra-low entry height allows seamless wheelchair transfers, and upholstery is MRSA-resistant, antibacterial, fire-retardant, and available in 14 standard colours.

We are confident in the materials and build quality of our products and as such offer a comprehensive parts and labour warranty to cover manufacturing or material defects of 7 years for the frame and 3 years for moving parts.

Code	Mechanism	SWL	Surface	Height Adj	Backrest Raises/Lowers
BO01E	Electric	320kg	190 x 102 cm	50-95 cm	-20° to +85°



Accessories	Code
Footswitch	FO01
Extra Wide Upholstery- large format couches only (120cm)	EW3
Breathing Hole	BR02
Detachable Head Rest	DHR03
Paper Roll Holder	PR01
Heavy Duty Adjustable Arm Rests	BE01
Battery Operating System	BOS
Cotsides (painted, chrome or padded)	




Find our products on the NHS Supply Chain Catalogue



14 Standard Colours



BEIGE



CITRUS GREEN



APPLE MINT



CANDY




RASPBERRY CRUSH



DAMSON



PILLAR BOX RED




WHITE



ASH GREY



GUN METAL GREY




BLACK



DARK BLUE



MID BLUE



SKY BLUE