

Physio Plinths



PHYS+02E/H, PHYS+03E/H, PHYS+04E/H, PHYS+05E/H

The Physio+ range of rehabilitation plinths is precision-engineered to support treatment scenarios where stability, adjustability, and clinician access are essential.

Each version ranging from the 2-section to the 5-section models offers identical performance credentials. All feature gas-strut assisted backrests, retractable easy-clean castors for seamless mobility and secure positioning for mobile hoist access. Every unit includes a stitched ergonomic breathing hole, a paper roll holder, adjustable leveling foot for uneven surfaces, and an extra-low entry height for wheelchair accessibility.


The 5-section model enhances treatment flexibility with added clinician access through adjustable head, arm supports, and a long breathing hole especially beneficial for face-down therapies and spinal mobilisation.

The 2-section version offers a streamlined, space-saving format that still delivers clinician access and robust support, ideal for standard examination or therapy workflows.

Code	Mechanism	SWL	Surface	Footrest Adj	Backrest Adj
PHYS+02E	Electric	260kg	186 x 63 cm	N/a	-20° to 80°
PHYS+02H	Hydraulic	225kg	186 x 63 cm	N/a	-20° to 80°
PHYS+03E	Electric	260kg	186 x 63 cm	80°	-20° to 80°
PHYS+03H	Hydraulic	225kg	186 x 63 cm	80°	-20° to 80°
PHYS+04E	Electric	260kg	186 x 63 cm	80°	-20° to 80°
PHYS+04H	Hydraulic	225kg	186 x 63 cm	80°	-20° to 80°
PHYS+05E	Electric	260kg	186 x 63 cm	80°	-20° to 80°
PHYS+05H	Hydraulic	225kg	186 x 63 cm	80°	-20° to 80°



Accessories	Code
Handset	HND00L
Extra Wide Upholstery (75cm)	EW1
Foot Protection Cover	FPO1
Cotsides (BA02/3 Only painted, chrome or padded)	
Positioning Wedge Set	WEDGESET



Find our products on the NHS Supply Chain Catalogue

14 Standard Colours



Physio Plinths




PHYS+BOE, PHYS+MTE

Engineered for rehabilitation settings such as physiotherapy, neurology, and head-injury care, the Physio+ Two-Section BO Plinth offers exceptional clinical access to support effective treatment. The plinth includes a stitched keyhole breathing aperture with bung, paper-roll holder, footswitch and a heavy-duty adjustable foot to ensure stability on uneven floors.

The Mat Table variant of the Bobath rehabilitation plinth offers the same robust performance and platform stability, including the 320kg electric capacity, Linak-driven lift, and high-grade hygienic upholstery. It similarly features gas-assisted positioning sections, footswitch control, paper-roll holder, retractable castors, and adjustable feet, it's ideal for extended physio routines or patient assessments.

Both units deliver exceptional durability and positioning flexibility for physiotherapy and rehabilitation environments. They offer enriched clinician access and patient comfort, all while meeting stringent infection-control standards and delivering long-term value through robust UK engineering.

Code	Mechanism	SWL	Surface	Footrest Adj	Backrest Adj
PHYS+BOE	Electric	320kg	186 x 102 cm	N/a	-20° to 80°
PHYS+MTE	Electric	320kg	186 x 102 cm	N/a	N/a




Accessories	Code
Handset	HND00L
Extra Wide Upholstery- large format couches only (120cm)	EW3
Cotsides (painted, chrome or padded)	
Positioning Wedge Set	WEDGESET



14 Standard Colours




BEIGE



CITRUS GREEN




APPLE MINT




CANDY




RASPBERRY CRUSH



DAMSON



PILLAR BOX RED




WHITE



ASH GREY



GUN METAL GREY



BLACK



DARK BLUE



MID BLUE



SKY BLUE