

Practice Plinths



PP02

The Medi-Plinth Practice Plinth (fixed-height 2-section couch) presents outstanding value without compromising durability or functionality, making it a preferred choice across varied medical settings. Constructed with a robust epoxy powder-coated steel frame, it delivers dependable performance in primary care, physiotherapy, and general clinical rooms. A standard paper roll holder, which can be positioned for clinician convenience, further enhances its practicality.

Upholstery is designed with clinic safety and hygiene in mind: seamless-film gussets are used between sections to minimise dirt traps, and the material is MRSA-resistant, fire-retardant, and available in 14 contemporary vinyl colours to suit any given interior scheme.

The Practice Plinth is backed by a three-year structural warranty (excluding upholstery), reflecting long-term product confidence and commitment to aftercare. This Plinth is delivered flat-pack for easier transport and access into restricted areas.



| Code | Mechanism | SWL | Surface | Height Adj | Raise/Lower Sections |
|------|--------------|-------|-------------|------------|----------------------|
| PP02 | Fixed Height | 180kg | 190 x 63 cm | n/a | 55° |

