

IV Poles



IVDST001/002/003/004, DST003

Medi-Plinth's Shuttleworth IV Poles are CE marked Class I Medical Devices, engineered for safe and reliable fluid administration in a range of healthcare settings. Built for stability and ease of use, these drip stands feature a quick-clip drop catcher, two or four hook configurations, and handle on models with a weighted base for enhanced manoeuvrability.

To reduce transport costs and environmental impact, IV poles are supplied flat-packed, while models with bases include five ESD swivel castors, two of which lock for added control.

A two-hook wall-mounted option is also available for space-limited areas.

Fully compliant with BS ISO 15375:2010, the Shuttleworth IV Pole range is a dependable, flexible solution for modern clinical care.



Code	Base Diameter	Adjustable Height	Hooks
IVDST001	Ø 71 cm	129-204cm	2
IVDST002	Ø 71 cm	129-204cm	2
IVDST003	Ø 71 cm	129-204cm	4
IVDST004	Ø 71 cm	129-204cm	4
DST003	N/a	60.8-120cm	2

Configuration

IVDST001: Two Hook IV Drip Stand, Stainless Steel, Flat Packable

IVDST002: Two Hook IV Drip Stand, Stainless Steel with Weighted Base, Flat Pack

IVDST003: Four Hook IV Drip Stand, Stainless Steel, Flat Packable

IVDST004: Four Hook IV Drip Stand, Stainless Steel with Weighted Base, Flat Pack

DST003: Two Hook Drip Stand, Stainless Steel, Wall Mounted