

Paediatric Cots



CHILD ASSURE Paediatric Cots

The Child Assure series offers purpose-built paediatric cots that prioritise child comfort, caregiver access, and clinical safety in busy neonatal and paediatric care settings.

The standard Child Assure Cot features a fixed height of 740 mm, an epoxy-coated steel frame with corner bumpers, and transparent side rails that offer 360° visibility while meeting strict entrapment safety standards. It offers a Trendelenburg tilt of $\pm 13^\circ$ via a hydraulic lever.

The Child Assure Care Cot extends functionality with electric height adjustability (62–92 cm) and electric Trendelenburg/reverse Trendelenburg tilt up to $\pm 16^\circ$. It includes clinician-operated handset controls, a quick-release backrest for emergency lowering, sturdy side rails with seven-position sliding adjustment.

For intensive care needs, the Child Assure Intensive Cot features four independent compact HPL panels for flexible patient positioning, hydraulic tilt functionality ($\pm 7^\circ / \pm 18^\circ$), and added hygiene and durability components such as an ABS-covered base with oxygen bottle storage.

Code	Dims	Back Section Dims	Fixed Section Dims	Leg Section Dims	SWL
CHILDASSURE-001	86 x 136.5	N/a	N/a	N/a	85kg
CHILDASSURE-002	157.5 x 155cm	N/a	N/a	N/a	85kg
CHILDASSURE-003	85.5 x 154cm	66 x 55cm	70.5 x 9.5cm	66 x 25cm	90kg

Code	Foot Section Dims	Backrest Tilt	Leg rest Angle	Height Adj	Trendelenburg
CHILDASSURE-001	N/a	N/a	N/a	Fixed 75cm	$\pm 13^\circ$
CHILDASSURE-002	N/a	N/a	N/a	58-90cm	$\pm 16^\circ$
CHILDASSURE-003	55 x 32.5cm	72°	36°	56-90cm	$\pm 7^\circ / \pm 18^\circ$

Configuration

CHILDASSURE-001: Child Assure (Paediatric Cot) - Standard Configuration

CHILDASSURE-002: Child Assure-Care (Paediatric Cot) - Standard Configuration

CHILDASSURE-003: Child Assure-Intensive (Paediatric Cot) - Standard Configuration

Accessories	Code
IV Pole for Paediatric Cot	CHILD-AC001
Polyurethane Mattress with Cover (Child Assure)	CHILD-AC002
Polyurethane Mattress with Cover (Child Assure-Care/Intensive)	CHILD-AC003

All models include integrated accessories and safety-oriented design, such as IV and mosquito net sockets, built-in bumpers, and sealed surfaces, tailored to infection control and child safety. The range meets Regulation (EU) 2017/745 compliance and is built to serve the unique demands of pediatric environments.

