

Easy Clean Overbed Tables



ECOBT01/02/03/04/05/06/07

Medi-Plinth Easy Clean Overbed Tables are designed to deliver dependable, hygienic, and versatile surface space for use in hospitals, care homes, and clinical facilities. Manufactured to BS 2483:1977 standards, they combine robust steel frames with a high-quality top surface that is easy to sanitise and withstands regular cleaning with NHS-approved solutions.

Available in white or beech finishes and with flat, lipped, or recessed top designs, they accommodate varied clinical and patient needs while maintaining a professional, modern appearance. With a safe working load of 25 kg, anti-entrapment safety features, and FIRA certification, Medi-Plinth overbed tables deliver dependable performance and patient-focused design backed by a 3-year manufacturer warranty.

The ECOBT07 Tilting Table offers excellent stability while maintaining a minimal footprint, making it ideal for use in confined clinical spaces. The integrated anti-slip retainer ensures items remain secure during use, while 360° swivel castors allow effortless positioning and manoeuvrability. Featuring a smooth height adjustment and a tilting tabletop of up to $\pm 33^\circ$.

Code	Total Width	Adjustable Height	Total Depth
ECOBT01+COLOUR CODE	90cm	79-118.5cm	40cm
ECOBT02+COLOUR CODE	90cm	79-118.5cm	40cm
ECOBT03+COLOUR CODE	90cm	79-118.5cm	40cm
ECOBT04+COLOUR CODE	90cm	79-118.5cm	40cm
ECOBT05+COLOUR CODE	90cm	79-118.5cm	40cm
ECOBT06+COLOUR CODE	90cm	79-118.5cm	40cm
ECOBT07	68cm	70-103.6cm	47.7cm

Configuration

- ECOBT01: 'C' leg base with flat top
- ECOBT02: 'C' leg base with upstand top
- ECOBT03: 'C' leg base with recessed top
- ECOBT04: Straight leg base with flat top
- ECOBT05: Straight leg base with upstand top
- ECOBT06: Straight leg base with recessed top
- ECOBT07: 'H' Shaped Base, Tilting Top - White



-LO



-GMG



-W