



PLINTHS & TREATMENT CHAIRS

Emphasizing Patient Comfort & Clinician Care

CONTENTS

» Ensignia Range	3
» 2 & 3 Section Plinths	5
» Changing Tables	7
» Drop End Plinths	9
» Multi-Purpose Plinths	11
» Phlebotomy Chair	13
» Treatment/ Plaster Chair	15
» Podiatry Chair	17
» Gynaecology	19
» Extra Width Plinths	23
» Echocardiography	27
» Outpatient Plinths	29
» Day Care Couches	31
» Medical Chairs	33
» Rehabilitation	37
» Imaging	47
» Fixed Height Plinths	51
UPHOLSTERY Colour Chart	53

Ensignia Range

E2/E3



Structure

Epoxy coated steel frame.
Industry-leading Linak 'High Reliability' actuators.
Electric lifting system provides optimum stability and supports weights of up to 260 kg.
Electric lifting with handset controls.
Easy clean, 100mm braked.
transportation castors provide excellent manoeuvrability.
Paper roll holder fitted as standard.



Easy Cleaning

Upholstery compatible with all NHS approved cleaning solutions such as Chlor-Clean, Haz-Tab, Milton and Tristel, plus Clinell wipes and other alcohol based wipes.



Comfort

Soft, durable vinyl to ensure patient comfort.
Patient-centric design enables enhanced positioning for improved comfort during examinations.



Safety

Anti-entrapment Suspa gas strut operated backrest/ Footrest. MRSA resistant anti-bacterial & fire retardant upholstery.

Comprehensive parts and labour warranty covering manufacturing or material defects of seven years for the frame and three years for moving parts.

Complies with all applicable medical regulations.



Intended Use

Engineered with precision, these plinths offer exceptional versatility, enabling a wide range of treatments, examinations, and therapies.

Key features include smooth electric height adjustment, adjustable sections for optimal patient positioning, and durable, easy-to-clean upholstery for maximum hygiene. The robust yet stylish design is complemented by transport castors, ensuring effortless mobility in busy clinical settings.

Specifications

ENSIGNIA RANGE

ENSIGNIA: TWO SECTION

Ensignia Two Section Plinth: **E2**

Patient surface	186 x 63 cm
Height adjustment	46 cm - 91cm
Backrest raises/lowers	-20° to +85°
Electric max. load	260 kg

Electric Features:

Power supply	110-240v / 50-60 Hz
Maximum current intake	5A, 100-240V
Protection indicator	IPX4
Protection class	Class 2.B
Electric shock protection	TYPE B

ENSIGNIA: THREE SECTION

Ensignia Three Section Plinth: **E3**

Patient surface	186 x 63 cm
Height adjustment	46 cm - 91cm
Backrest raises/lowers	-20° to +85°
Electric max. load	260 kg
Footrest raises/lowers	-20° to +85°

Optional Accessories:

- « Footswitch
- « Extra Wide Upholstery (75 cm)
- « Breathing Hole
- « Foot Protection Cover
- « Detachable Head Rest
- « Cotsides
- « Clean Edge Clover



			
BEIGE	CITRUS GREEN	APPLE MINT	CANDY
			
WHITE	ASH GREY	GUN METAL GREY	BLACK
			
RASPBERRY CRUSH	DAMSON	PILLAR BOX RED	DARK BLUE
			
MID BLUE	SKY BLUE		

2 & 3 Section Plinths

BA02E/H, BA03E/H



Structure

Epoxy coated steel frame.
Industry-leading Linak 'High Reliability' actuators.
Electric lifting system provides optimum stability and supports weights of up to 260 kg.
Hydraulic lifting system provides optimum stability and supports weights of up to 225kg.
Electric lifting with handset controls.
Electric backrest, foot and/ or tilting seat available.
Paper roll holder fitted as standard.



Easy Cleaning

Upholstery compatible with all NHS approved cleaning solutions such as Chlor-Clean, Haz-Tab, Milton and Tristel, plus Clinell wipes and other alcohol based wipes.



Comfort

Soft, durable vinyl to ensure patient comfort.
Patient-centric design enables enhanced positioning for improved comfort during examinations.



Safety

Anti-entrapment Suspa gas strut operated backrest/footrests. MRSA resistant anti-bacterial & fire retardant upholstery.

Comprehensive parts and labour warranty covering manufacturing or material defects of seven years for the frame and three years for moving parts.



Intended Use

Designed to provide exceptional functionality and comfort for both patients and clinicians. Its adjustable sections and a robust construction, mean these couches support a wide range of medical procedures, from general examinations to specialist treatments.

Key features include smooth height adjustments, durable upholstery for easy cleaning, and ergonomic design for patient comfort and clinician convenience.

Specifications

2 & 3 SECTION PLINTHS

TWO SECTION

2 Section Electric Plinth: **BA02E**

2 Section Hydraulic Plinth: **BA02H**

Patient surface	186 x 63 cm
Height adjustment	46cm - 91cm
Backrest raises/lowers	-20° to +85°
Electric max. load	260kg
Hydraulic max. load	225kg

Electric Features:

Power supply	110-240v / 50-60
Maximum current intake	5A, 100-240V
Protection indicator	IPX4
Protection class	Class 2.B
Electric shock protection	TYPE B



THREE SECTION

3 Section Electric Plinth: **BA03E**

3 Section Hydraulic Plinth: **BA03H**

Patient surface	186 x 63 cm
Height adjustment	46cm - 91cm
Backrest raises/lowers	-20° to +85°
Electric max. load	260 kg
Hydraulic max. load	225 kg
Footrest raises/lowers	-20 - 85°
Trendelenburg (Tilting Versions)	-15°

Optional Accessories:

- « Footswitch
- « Extra Wide Upholstery (75 cm)
- « Breathing Hole
- « Foot Protection Cover
- « Detachable Head Rest
- « Cotsides
- « Tilting Seat (3 Sections only)
- « Electric Foot (3 Sections only)
- « Electric Backrest
- « Clean Edge Cover



Changing Tables

BA01E/BA01H



Structure

Epoxy coated steel frame.
Industry-leading Linak 'High Reliability' actuators.
Electric lifting system provides optimum stability and supports weights of up to 260 kg.
Hydraulic lifting system provides optimum stability and supports weights of up to 225kg.
Electric lifting with handset controls.
Available in sizes 1000mm, 1200mm, 1500mm or 1800mm.



Easy Cleaning

Upholstery compatible with all NHS approved cleaning solutions such as Chlor-Clean, Haz-Tab, Milton and Tristel, plus Clinell wipes and other alcohol based wipes.



Comfort

Soft, durable vinyl to ensure patient comfort.
Patient-centric design enables enhanced positioning for improved comfort.



Safety

MRSA resistant anti-bacterial & fire retardant upholstery.

Comprehensive parts and labour warranty covering manufacturing or material defects of seven years for the frame and three years for moving parts.

Complies with all applicable medical regulations.



Intended Use

Designed with durability and practicality in mind, these tables are ideal for use in medical facilities, schools, and care environments.

Featuring smooth height adjustment for ease of use, sturdy construction for enhanced safety, and easy-to-clean upholstery, our changing tables are suitable for both children and adults requiring assisted care.



Please note cotsides are not included as standard.

Specifications

CHANGING TABLES

PATIENT SURFACE OPTIONS

BA01E1000	1000 x 630 mm
BA01E1200	1200 x 630 mm
BA01E1500	1500 x 630 mm
BA01E1800	1800 x 630 mm

Electric Features:

Power supply	110-240v / 50-60
Maximum current intake	5A, 100-240V
Protection indicator	IPX4
Protection class	Class 2.B
Electric shock protection	TYPE B

ELECTRIC/ HYDRAULIC MECHANISM

Changing Table Electric: **BA01E**

Changing Table Hydraulic: **BA01H**


Height adjustment	46 cm - 91cm
Electric max. load	260 kg
Hydraulic max. load	225 kg

Optional Accessories:

- « Cotsides
- « Paper roll holder
- « Footswitch
- « Extra-wide upholstery (750mm)
- « Clean Edge Cover

COTSIDE OPTIONS:

- « Cotside Chromed Plated - Pair: **CO01C-P**
- « Cotside Chrome Plated (Single): **CO01C-S**
- « Cotside Padded Loop Chrome Plated - Pair: **CP01C-P**
- « Cotside Padded Loop Chrome Plated (Single): **CP01C-S**
- « Cotside Painted Padded- Pair: **CP01-P**
- « Cotside Painted Padded (Single): **CP01-S**
- « Cotside Painted- Pair: **CO01-P**
- « Cotside Painted (Single): **CO01-S**



Please note cotsides are not included as standard.

BEIGE	CITRUS GREEN	APPLE MINT	CANDY	RASPBERRY CRUSH	DAMSON	PILLAR BOX RED
WHITE	ASH GREY	GUN METAL GREY	BLACK	DARK BLUE	MID BLUE	SKY BLUE

Drop End Plinths

DE03E/DE03H



Structure

Epoxy coated steel frame.
Industry-leading Linak 'High Reliability' actuators.
Electric lifting system provides optimum stability and supports weights of up to 260 kg.
Hydraulic lifting system provides optimum stability and supports weights of up to 225kg.
Electric lifting with handset controls.
Optional electric back and/or tilting seat.



Easy Cleaning

Upholstery compatible with all NHS approved cleaning solutions such as Chlor-Clean, Haz-Tab, Milton and Tristel, plus Clinell wipes and other alcohol based wipes.



Comfort

Soft, durable vinyl to ensure patient comfort.
Patient-centric design enables enhanced positioning for improved comfort during examinations.



Safety

Anti-entrapment Suspa gas strut operated back/footrest. MRSA resistant anti-bacterial & fire retardant upholstery.

Comprehensive parts and labour warranty covering manufacturing or material defects of seven years for the frame and three years for moving parts.

Complies with all applicable medical regulations.



Intended Use

The Medi-Plinth Drop End Plinth is a practical and versatile solution designed to support a wide range of examinations and treatments. Featuring an adjustable drop-end section, this plinth allows for enhanced patient positioning and improved accessibility, making it ideal for procedures requiring seated or reclined postures.

Specifications

DROP END PLINTHS

DROP END PLINTH

3 Section Electric Drop End Plinth: **DE03E**

Patient surface	186 x 63 cm
Height adjustment	46 cm - 91 cm
Backrest raises/lowers	-20 to +85°
Trendelenburg (tilting versions)	-15°
Electric max. load	260 kg
Footrest lowers	-85°

Electric Features:

Power supply	110-240v / 50-60
Maximum current intake	5A, 100-240V
Protection indicator	IPX4
Protection class	Class 2.B
Electric shock protection	TYPE B



DROP END PLINTH

3 Section Hydraulic Drop End Plinth: **DE03H**

Patient surface	186 x 63 cm
Height adjustment	46 cm - 91 cm
Backrest raises/lowers	-20 to + 85°
Hydraulic max. load	225 kg
Footrest lowers	-85°

Optional Accessories:

- « Footswitch
- « Extra Wide Upholstery (750mm)
- « Breathing Hole
- « Foot Protection Cover
- « Detachable Head Rest
- « Adjustable Arm Rests
- « Tilting Seat (3 Sections only)
- « Electric Backrest (3 Sections only)
- « Stainless Steel Accessory Rails, for the use of gynae or arm support attachments
- « Clean Edge Cover



Multi-Purpose Plinths

MP03E/MP03H



Structure

Epoxy coated steel frame.
Industry-leading Linak 'High Reliability' actuators.
Electric lifting system provides optimum stability and supports weights of up to 260 kg.
Hydraulic lifting system provides optimum stability and supports weights of up to 225kg.
Electric lifting with handset controls.
Includes: 1x Phlebotomy arm, 2x Gynaecology stirrups & 1x Quick release plaster treatment table.
Optional electric back and/or tilting seat.



Easy Cleaning

Upholstery compatible with all NHS approved cleaning solutions such as Chlor-Clean, Haz-Tab, Milton and Tristel, plus Clinell wipes and other alcohol based wipes.



Comfort

Soft, durable vinyl to ensure patient comfort.
Patient-centric design enables enhanced positioning for improved comfort during examinations.



Safety

Anti-entrapment Suspa gas strut operated back/footrest. MRSA resistant anti-bacterial & fire retardant upholstery.

Comprehensive parts and labour warranty covering manufacturing or material defects of seven years for the frame and three years for moving parts.

Complies with all applicable medical regulations.



Intended Use

The Multi-Purpose Plinth is a flexible and reliable solution designed to meet the needs of a wide range of medical procedures. Whether for general examinations or minor treatments, this plinth delivers exceptional performance in any clinical setting.

With smooth height adjustment, adjustable sections, and high-quality, easy-to-clean upholstery, the Multi-Purpose Plinth ensures maximum comfort for patients and convenience for clinicians.

Specifications

MULTI-PURPOSE PLINTHS

MULTI-PURPOSE PLINTH

3 Section Electric Multi-Purpose Plinth: **MP03E**

Patient surface	177 x 63 cm
Height adjustment	46 cm - 91 cm
Backrest raises/lowers	-20° to +85°
Trendelenburg (tilting versions)	-15°
Electric max. load	260 kg
Footrest lowers	-85°

Electric Features:

Power supply	110-240v / 50-60
Maximum current intake	5A, 100-240V
Protection indicator	IPX4
Protection class	Class 2.B
Electric shock protection	TYPE B

MULTI-PURPOSE PLINTH

3 Section Hydraulic Multi-Purpose Plinth: **MP03H**

Patient surface	177 x 63 cm
Height adjustment	46 cm - 91 cm
Backrest raises/lowers	-20° to +85°
Hydraulic max. load	225 kg
Footrest lowers	-85°

Optional Accessories:

- « Footswitch
- « Extra Wide Upholstery (750mm)
- « Breathing Hole
- « Foot Protection Cover
- « Detachable Head Rest
- « Tilting Seat
- « Electric Backrest
- « Paper Roll Holder
- « Clean Edge Cover



			
BEIGE	CITRUS GREEN	APPLE MINT	CANDY
			
WHITE	ASH GREY	GUN METAL GREY	BLACK
			
RASPBERRY CRUSH	DAMSON	PILLAR BOX RED	DARK BLUE
			
MID BLUE	SKY BLUE		

Phlebotomy

PH01E/H, PH02E/H



Structure

Epoxy coated steel frame.
Industry-leading Linak 'High Reliability' actuators.
Electric lifting system provides optimum stability and supports weights of up to 260 kg.
Hydraulic lifting system provides optimum stability and supports weights of up to 225kg.
Electric lifting with handset controls.
Two phlebotomy arms as standard.
Available in split leg or single leg models.
Optional electric back and/or tilting seat.



Easy Cleaning

Upholstery compatible with all NHS approved cleaning solutions such as Chlor-Clean, Haz-Tab, Milton and Tristel, plus Clinell wipes and other alcohol based wipes.



Comfort

Soft, durable vinyl to ensure patient comfort.
Patient-centric design enables enhanced positioning for improved comfort during examinations. Single or split-leg versions available.



Safety

Anti-entrapment Suspa gas strut operated back/footrest. MRSA resistant anti-bacterial & fire retardant upholstery.
Comprehensive parts and labour warranty covering manufacturing or material defects of seven years for the frame and three years for moving parts.

Complies with all applicable medical regulations.



Intended Use

The Phlebotomy chair is designed to provide optimal comfort and support during blood collection procedures. Featuring a stable and adjustable design, this plinth allows for precise positioning of patients, ensuring easy access for phlebotomists while maintaining patient comfort throughout the procedure.

With smooth height adjustments, durable, easy-to-clean upholstery, and a strong, reliable frame, the Plinth is perfect for clinics, blood donation centers, and medical offices. Its compact design makes it ideal for high-traffic environments, delivering efficient performance without compromising on patient care.

Specifications

PHLEBOTOMY CHAIRS

PHLEBOTOMY CHAIR

Phlebotomy Chair Hydraulic Single Leg: **PH02H**

Phlebotomy Chair Hydraulic Split Leg: **PH01H**

Patient surface	186 x 63 cm
Height adjustment	46 cm - 91 cm
Backrest raises/lowers	-20° to +85°
Hydraulic max. load	225 kg
Footrest lowers	- 85°

Electric Features:

Power supply	110-240v / 50-60
Maximum current intake	5A, 100-240V
Protection indicator	IPX4
Protection class	Class 2.B
Electric shock protection	TYPE B



PHLEBOTOMY CHAIR

Phlebotomy Chair Electric Single Leg: **PH02E**

Phlebotomy Chair Electric Split Leg: **PH01E**

Patient surface	186 x 63 cm
Height adjustment	46 cm - 91 cm
Backrest raises/lowers	-20° to +85°
Trendelenburg (tilting versions)	-15°
Electric max. load	260 kg
Footrest lowers	- 85°

Optional Accessories:

- « Footswitch
- « Extra Wide Upholstery (750mm)
- « Breathing Hole
- « Foot Protection Cover
- « Detachable Head Rest
- « Tilting Seat (3 Sections only)
- « Electric Backrest
- « Paper Roll Holder
- « Clean Edge Cover



Treatment/Plaster

PL01E/H, PL02E/H



Structure

Epoxy coated steel frame.
Industry-leading Linak 'High Reliability' actuators.
Electric lifting system provides optimum stability and supports weights of up to 260 kg.
Hydraulic lifting system provides optimum stability and supports weights of up to 225kg.
Electric lifting with handset controls.
Optional electric back and/or tilting seat.
2x Quick release treatment tables



Easy Cleaning

Upholstery compatible with all NHS approved cleaning solutions such as Chlor-Clean, Haz-Tab, Milton and Tristel, plus Clinell wipes and other alcohol based wipes.



Comfort

Soft, durable vinyl to ensure patient comfort.
Patient-centric design enables enhanced positioning for improved comfort during examinations. Single or split-leg versions available.



Safety

Anti-entrapment Suspa gas strut operated back/footrest. MRSA resistant anti-bacterial & fire retardant upholstery.
Comprehensive parts and labour warranty covering manufacturing or material defects of seven years for the frame and three years for moving parts.

Complies with all applicable medical regulations.



Intended Use

The Medi-Plinth Treatment/Plaster Chair is engineered to provide exceptional support and stability during a wide range of treatments, including plaster casting, and wound care. Its strong frame and durable construction ensure it can withstand the demands of everyday use in busy healthcare settings.

Featuring adjustable height, a smooth and easy-to-clean surface, and plaster arm tables, this plinth ensures both patient comfort and clinician convenience.

Specifications

TREATMENT/PLASTER CHAIRS

TREATMENT/PLASTER CHAIR

Treatment/Plaster Chair Electric - Split leg: **PL01E**

Treatment/Plaster Chair Hydraulic - Split leg: **PL01H**

Patient surface	186 x 63 cm
Height adjustment	46 cm - 91cm
Backrest raises/lowers	-20° to +85°
Hydraulic max. load	225 kg
Footrest lowers	-85°
Trendelenburg (tilting versions)	-15°

Electric Features:

Power supply	110-240v / 50-60
Maximum current intake	5A, 100-240V
Protection indicator	IPX4
Protection class	Class 2.B
Electric shock protection	TYPE B



TREATMENT/PLASTER CHAIR

Treatment/Plaster Chair Hydraulic - Single leg: **PL02H**

Treatment/Plaster Chair Electric - Single leg: **PL02E**

Patient surface	186 x 63 cm
Height adjustment	46 cm - 91cm
Backrest raises/lowers	-20° to +85°
Trendelenburg (tilting versions)	-15°
Electric max. load	260 kg
Footrest lowers	-85°

Optional Accessories:

- « Footswitch
- « Extra Wide Upholstery (75 cm)
- « Breathing Hole
- « Foot Protection Cover
- « Detachable Head Rest
- « Tilting Seat
- « Electric Backrest
- « Paper Roll Holder
- « Clean Edge Cover



Podiatry

BAC06E/BTC06E

**Structure**

Epoxy coated steel frame.
Industry-leading Linak 'High Reliability' actuators.
Electric lifting system provides optimum stability and supports weights of up to 260 kg.
Electric lifting with handset controls as standard.
Two fully adjustable armrests.
Optional electric back and/or tilting seat.

**Easy Cleaning**

Upholstery compatible with all NHS approved cleaning solutions such as Chlor-Clean, Haz-Tab, Milton and Tristel, plus Clinell wipes and other alcohol based wipes.

**Comfort**

Soft, durable vinyl to ensure patient comfort.
Patient-centric design enables enhanced positioning for improved comfort during examinations.

**Safety**

Anti-entrapment Suspa gas strut operated back/footrest. MRSA resistant anti-bacterial & fire retardant upholstery.
Comprehensive parts and labour warranty covering manufacturing or material defects of seven years for the frame and three years for moving parts.

Complies with all applicable medical regulations.

**Intended Use**

Medi-Plinth's Podiatry Chairs are designed to offer superior comfort, stability, and adjustability for podiatric examinations and treatments. Engineered for both patient comfort and clinician efficiency, these chairs feature fully adjustable positions, allowing for optimal positioning during foot care procedures.

Their practical design ensures ease of access for clinicians while enhancing patient comfort during extended treatments.

Specifications

PODIATRY CHAIRS

PODIATRY CHAIR

Podiatry Chair Electric - Split leg: **BAC06E**

Patient surface	186 x 63 cm
Height adjustment	46 cm - 91 cm
Backrest raises/lowers	-20 - 85°
Footrest lowers	-85°
Electric max. load	260 kg

Electric Features:

Power supply	110-240v / 50-60
Maximum current intake	5A, 100-240V
Protection indicator	IPX4
Protection class	Class 2.B

PODIATRY CHAIR

Podiatry Chair Electric - Split leg, Tilting: **BTC06E**

Patient surface	186 x 63 cm
Height adjustment	46 cm - 91 cm
Backrest raises/lowers	-20 - 85°
Trendelenburg (tilting versions)	-15°
Electric max. load	260 kg
Footrest lowers	-85°

Optional Accessories:

- « Footswitch
- « Extra Wide Upholstery (75 cm)
- « Foot Protection Cover
- « Detachable Head Rest
- « Tilting Seat
- « Electric Backrest
- « Paper Roll Holder
- « Clean Edge Cover



BEIGE	CITRUS GREEN	APPLE MINT	CANDY
WHITE	ASH GREY	GUN METAL GREY	BLACK
RASPBERRY CRUSH	DAMSON	PILLAR BOX RED	DARK BLUE
MID BLUE	SKY BLUE		

Gynaecology

GY02ET/GY02ES/GY02EF



Structure

Epoxy coated steel frame.
Industry-leading Linak 'High Reliability' actuators.
Electric lifting system provides optimum stability and supports weights of up to 260 kg.
Electric lifting with handset controls.
Removable leg section enables couch to be used as a standard two section plinth.
The electric seat enables Trendelenburg positioning.



Easy Cleaning

Upholstery compatible with all NHS approved cleaning solutions such as Chlor-Clean, Haz-Tab, Milton and Tristel, plus Clinell wipes and other alcohol based wipes.



Comfort

Soft, durable vinyl to ensure patient comfort. Patient-centric design enables enhanced positioning for improved comfort during examinations. Optional accessory of either troughs, footpads or stirrups included.



Safety

Anti-entrapment Suspa gas strut operated backrest. MRSA resistant anti-bacterial & fire retardant upholstery.
Comprehensive parts and labour warranty covering manufacturing or material defects of seven years for the frame and three years for moving parts.
Complies with all applicable medical regulations.



Intended Use

Combining comfort, stability, and functionality, this plinth provides a secure and supportive surface for patients during examinations, treatments, and minor procedures.

Featuring adjustable height, a comfortable and easy-to-clean surface, and robust construction, the GP Gynae Plinth ensures optimal positioning for both patients and clinicians. Its versatile design makes it ideal for a wide range of medical applications, including gynaecological, obstetric, and general examinations.

Specifications

GP GYNAE PLINTH

GP GYNAE: THREE SECTION

With Foot pads: **GY02EF**

With Stirrups: **GY02ES**

With Troughs: **GY02ET**

**Removeable leg section included as standard.
Can be used as a standard two section plinth.**



Patient surface	173 x 63 cm
Height adjustment	50–95 cm
Backrest raises/lowers	-20° to +85°
Trendelenburg	-15°
Electric max. load	260 kg

Power supply	110-240v / 50-60 Hz
Maximum current intake	5A, 100-240V
Protection indicator	IPX4
Protection class	Class 2.B
Electric shock protection	TYPE B

Optional Accessories:

- « Footswitch
- « Extra Wide Upholstery (75 cm)
- « Breathing Hole
- « Foot Protection Cover
- « Detachable Head Rest
- « Electric Backrest
- « Paper Roll Holder
- « Fold Down Arm Loops
- « Troughs, Stirrups or Foot Pads



Please note Fold Down Arm Loops are not included as standard.



						
BEIGE	CITRUS GREEN	APPLE MINT	CANDY	RASPBERRY CRUSH	DAMSON	PILLAR BOX RED
						
WHITE	ASH GREY	GUN METAL GREY	BLACK	DARK BLUE	MID BLUE	SKY BLUE

Gynaecology

GY04ET/GY04EF/GY04ES



Structure

Epoxy coated steel frame.
100% usable patient surface.
Features three electric motors for height, tilt and backrest control.
Industry-leading Linak 'High Reliability' actuators.
Stainless steel tray included as standard.
Electric lifting system provides optimum stability and supports weights of up to 260 kg.
Optional electric backrest available.
Optional accessory of either troughs, footpads or stirrups included.



Easy Cleaning

Upholstery compatible with all NHS approved cleaning solutions such as Chlor-Clean, Haz-Tab, Milton and Tristel, plus Clinell wipes and other alcohol based wipes.



Comfort

Soft, durable vinyl to ensure patient comfort.
Patient-centric design enables enhanced positioning for improved comfort during examinations.



Safety

Anti-entrapment Suspa gas strut operated backrest. MRSA resistant anti-bacterial & fire retardant upholstery.
Comprehensive parts and labour warranty covering manufacturing or material defects of seven years for the frame and three years for moving parts.
Complies with all applicable medical regulations.



Intended Use

The Medi-Plinth Ultra Gynae Plinth is designed to provide the highest level of comfort, support, and functionality for gynaecological examinations and procedures. Featuring an innovative design, this plinth offers adjustable positioning, ensuring optimal comfort for patients and ease of access for clinicians.

With smooth height adjustments, a comfortable, easy-to-clean surface, and strong, reliable construction, the Ultra Gynae Plinth is built for durability and long-lasting performance in busy clinical environments. for a variety of treatments.

Specifications

ULTRA GYNAE PLINTH

ULTRA GYNAE PLINTHS

With Foot Pads: **GY04EF**

With Stirrups: **GY04ES**

With Troughs: **GY02ET**



Patient surface	184 x 63 cm
Height adjustment	60 cm - 105.5 cm
Backrest raises	85°
Legrest raises	85°
Mid-section raises	85°
Trendelenburg	-20°
Electric max load	260 kg

Electric Features:

Power supply	110-240v / 50-60 Hz
Maximum current intake	5A, 100-240V
Protection indicator	IPX4
Protection class	Class 2.B
Electric shock protection	TYPE B

Optional Accessories:

- « Footswitch
- « Extra Wide Upholstery (75 cm)
- « Breathing Hole
- « Foot Protection Cover
- « Detachable Head Rest
- « Electric Backrest
- « Paper Roll Holder
- « Fold down arm loops
- « Troughs, Stirrups or Foot Pads



						
BEIGE	CITRUS GREEN	APPLE MINT	CANDY	RASPBERRY CRUSH	DAMSON	PILLAR BOX RED
						
WHITE	ASH GREY	GUN METAL GREY	BLACK	DARK BLUE	MID BLUE	SKY BLUE

Extra Wide Plinths

BO01E



Structure

Epoxy coated steel frame.
100% usable patient surface.
Excellent lithotomy access and foot space clinician.
Industry-leading Linak 'High Reliability' actuators.
Electric lifting system provides optimum stability and supports weights of up to 260 kg.
Electric lifting with handset controls.
Optional electric backrest available.
Mobile hoist clearance space beneath base.



Easy Cleaning

Upholstery compatible with all NHS approved cleaning solutions such as Chlor-Clean, Haz-Tab, Milton and Tristel, plus Clinell wipes and other alcohol based wipes.



Comfort

Soft, durable vinyl to ensure patient comfort. Its sturdy, wide top makes the couch safe for the therapist to treat the patient with ease and ensures patient comfort.



Safety

Anti-entrapment Suspa gas strut operated backrest. MRSA resistant anti-bacterial & fire retardant upholstery.
Comprehensive parts and labour warranty covering manufacturing or material defects of seven years for the frame and three years for moving parts.
Complies with all applicable medical regulations.



Intended Use

Designed to support neurological rehabilitation for individuals with conditions such as stroke, cerebral palsy, and other motor disorders. With adjustable sections and a sturdy frame, this plinth allows for versatile positioning, enhancing patient comfort and therapist control during treatment sessions.

Ideal for physiotherapists and rehabilitation centers, it offers an essential tool to facilitate motor recovery and improve functional outcomes.

Specifications

2 Section Bo-Plinth

2 Section Bo-Plinth

40" Bo-Plinth 2 Section Electric : **BO01E**

Electric Features:

Patient surface	190 x 102 cm
Height adjustment	50 cm - 95 cm
Backrest raises	85°
Backrest Lowers	-20°
Electric max load	320 kg

Power supply	110-240v / 50-60 Hz
Maximum current intake	5A, 100-240V
Protection indicator	IPX4
Protection class	Class 2.B
Electric shock protection	TYPE B
Backup lithium battery output	24V DC – 2.2 Ah

Optional Accessories:

- « Footswitch
- « Breathing Hole
- « Foot Protection Cover
- « Detachable Head Rest
- « Paper Roll Holder
- « Cotsides
- « Battery Operating System



Extra Wide Plinths

BARI0275/100, BARI0375/100, BARI0375DE/100DE



Structure

Available in 75 or 100cm wide.
Available as 2 or 3 section plinths.
Drop end versions also available.
Epoxy coated steel frame.
Industry-leading Linak 'High Reliability' actuators.
Electric lifting system provides optimum stability and supports weights of up to 320 kg.
Electric backrest as standard.



Easy Cleaning

Upholstery compatible with all NHS approved cleaning solutions such as Chlor-Clean, Haz-Tab, Milton and Tristel, plus Clinell wipes and other alcohol based wipes.



Comfort

Soft, durable vinyl to ensure patient comfort. Its sturdy, wide top makes the couch safe for the therapist to treat the patient with ease and ensures patient comfort.



Safety

MRSA resistant anti-bacterial & fire retardant upholstery.
Comprehensive parts and labour warranty covering manufacturing or material defects of seven years for the frame and three years for moving parts.
Complies with all applicable medical regulations.



Intended Use

These Bariatric Plinths are expertly designed to provide exceptional support and comfort for larger patients, while ensuring safety and ease of use for practitioners. Built with a sturdy, reinforced frame and a high weight capacity, these couches deliver durability and peace of mind in every healthcare setting.

Built to support a higher weight capacity without compromising on patient comfort, they are ideal for hospitals, rehabilitation centers, and healthcare environments focused on inclusive care.

Specifications

Bariatric Plinths

BARIATRIC: TWO SECTION

Bariatric 2 Section Plinth, 75cm: **BARI0275**

Bariatric 2 Section Plinth, 100cm: **BARI0100**

Patient surface (75cm)	186 x 75 cm
Patient surface (100cm)	186 x 100 cm
Height adjustment	50 cm - 95 cm
Backrest raises/lowers	-20° to +85°
Electric max. load	320 kg

Electric Features:

Power supply	110-240v / 50-60 Hz
Maximum current intake	5A, 100-240V
Protection indicator	IPX4
Protection class	Class 2.B
Electric shock protection	TYPE B

Optional Accessories:

- « Footswitch
- « Breathing Hole
- « Foot Protection Cover
- « Detachable Head Rest
- « Paper Roll Holder
- « Cotsides
- « Heavy Duty Adjustable Armrests

BARIATRIC: THREE SECTION

Bariatric 3 Section Plinth, 75cm: **BARI0375**

Bariatric 3 Section Plinth, 100cm: **BARI0100**

Bariatric 3 Section Plinth, 100cm, Drop End: **BARI0100DE**

Bariatric 3 Section Plinth, 75cm, Drop End: **BARI0375DE**

Patient surface (75cm)	186 x 75 cm
Patient surface (100cm)	186 x 100 cm
Height adjustment	50 cm - 95 cm
Backrest raises/lowers	-20° to +85°
Footrest raises/lowers	-20° to +85°
Footrest lowers (Drop end model only)	-85°
Electric max. load	320 kg



Echocardiography

ECG03E & ECG03T



Structure

Available with electric back rest and/or tilting seat. Industry-leading Linak 'High Reliability' actuators. Electric lifting system provides optimum stability and supports weights of up to 260 kg. Electric lifting with handset controls. Patient steadying handles fitted as standard. Extending operator seat fitted as standard.



Easy Cleaning

Upholstery compatible with all NHS approved cleaning solutions such as Chlor-Clean, Haz-Tab, Milton and Tristel, plus Clinell wipes and other alcohol based wipes.



Comfort

Soft, durable vinyl to ensure patient comfort. Its sturdy top makes the couch safe for the therapist to treat the patient with ease and ensures patient comfort.



Safety

Anti-entrapment Suspa gas strut operated back-rest/footrest. MRSA resistant anti-bacterial & fire retardant upholstery. Comprehensive parts and labour warranty covering manufacturing or material defects of seven years for the frame and three years for moving parts. Complies with all applicable medical regulations.



Intended Use

Designed with cardiology professionals in mind, the Medi-Plinth Echocardiography Plinth ensures optimal patient positioning and clinician access for precise diagnostics. Featuring an ergonomic cutaway design, this plinth provides effortless access to the thoracic region, enabling accurate and efficient echocardiographic examinations.

Specifications

Echocardiography Plinths

ECHOCARDIOGRAPHY

Echocardiography Plinth Electric -
Non-Tilting: **ECG03E**

Patient surface	205 x 74 cm
Height adjustment	46 cm - 91 cm
Backrest raises/lowers	-20° to +85°
Footrest lowers	-85°
Electric max. load	260 kg

ECHOCARDIOGRAPHY

Echocardiography Plinth Electric - Tilting:
ECG03T

Patient surface	205 x 74 cm
Height adjustment	46 cm - 91 cm
Backrest raises/lowers	-20° to +85°
Trendelenburg (tilting versions)	-15°
Electric max. load	260 kg
Footrest lowers	-85°



Optional Accessories:

- « Footswitch
- « Foot Protection Cover
- « Detachable Head Rest
- « Paper Roll Holder
- « Clean Edge Cover



Outpatient

OP02E/H, OP03E/H



Structure

Available with electric back rest and/ or electric foot.
Industry-leading Linak 'High Reliability' actuators. Electric lifting system provides optimum stability and supports weights of up to 260 kg. Hydraulic lifting system provides optimum stability and supports weights of up to 225kg. Fitted with transportation castors as standard. Fitted with chrome cotsides and BOS as standard.



Easy Cleaning

Upholstery compatible with all NHS approved cleaning solutions such as Chlor-Clean, Haz-Tab, Milton and Tristel, plus Clinell wipes and other alcohol based wipes.



Comfort

Soft, durable vinyl to ensure patient comfort. Its sturdy top makes the couch safe for the therapist to treat the patient with ease and ensures patient comfort.



Safety

Anti-entrapment Suspa gas strut operated back-rest/footrest as standard.
MRSA resistant anti-bacterial & fire retardant upholstery.
Comprehensive parts and labour warranty covering manufacturing or material defects of seven years for the frame and three years for moving parts.
Complies with all applicable medical regulations.



Intended Use

Medi-Plinth's 2 and 3 Section Outpatient Plinths are expertly designed to meet the demands of modern outpatient care, making them ideal for examinations, therapies, and minor procedures.

Key features include smooth height adjustments, durable upholstery, and multiple section configurations to enhance positioning. Equipped with high-quality transport castors, these plinths offer effortless mobility, ensuring flexibility and convenience in any clinical environment.

Specifications

2 & 3 Section Out Patient Plinths

TWO SECTION

Outpatients 2 Section Plinth, Hydraulic: **OP02H**

Outpatients 2 Section Plinth, Electric: **OP02E**

Patient surface	186 x 63 cm
Height adjustment	46 cm - 91cm
Backrest raises/ lowers	-20° to +85°
Electric max. load	260 kg
Hydraulic max. load	225 kg

Electric Features:

Power supply	110-240v / 50-60
Maximum current intake	5A, 100-240V
Protection indicator	IPX4
Protection class	Class 2.B
Electric shock protection	TYPE B



THREE SECTION

Outpatients 3 Section Plinth, Hydraulic: **OP03H**

Outpatients 3 Section Plinth, Electric: **OP03E**

Patient surface	186 x 63 cm
Height adjustment	46 cm - 91 cm
Backrest raises/ lowers	-20° to +85°
Electric max. load	260 kg
Hydraulic max. load	225 kg
Footrest raises/ lowers	-20° to +85°

Optional Accessories:

- « Footswitch (For Electric Versions)
- « Extra Wide Upholstery (75 cm)
- « Breathing Hole
- « Foot Protection Cover
- « Detachable Head Rest
- « Electric Foot (3 Sections Only)
- « Electric Backrest
- « IV Pole attachment



Day Care DC03E



Structure

Epoxy coated steel frame.
100% usable patient surface.
Industry-leading Linak 'High Reliability' actuators.
Electric lifting system provides optimum stability and supports weights of up to 260 kg.
Linked castor wheel base as standard.
Fold down sides/armrests as standard.



Intended Use

Medi-Plinth's Day Care Couches are designed to provide exceptional comfort and support for patients undergoing day treatments or short-stay procedures. These couches ensure a seamless experience for both patients and healthcare professionals.

With adjustable features, premium-quality upholstery, and robust construction, our day care couches are ideal for infusion therapies, recovery rooms, and general outpatient care.



Safety

Anti-entrapment Suspa gas strut operated back and footrest.
MRSA resistant anti-bacterial & fire retardant upholstery.
Comprehensive parts and labour warranty covering manufacturing or material defects of seven years for the frame and three years for moving parts.
Complies with all applicable medical regulations.



Easy Cleaning

Upholstery compatible with all NHS approved cleaning solutions such as Chlor-Clean, Haz-Tab, Milton and Tristel, plus Clinell wipes and other alcohol based wipes.



Comfort

Soft, durable vinyl to ensure patient comfort. Its sturdy, wide top makes the couch safe for the therapist to treat the patient with ease and ensures patient comfort.

Specifications

DAY CARE COUCH

DAY CARE

Day Care Couch Electric - Tilting: **DC03E**

Patient surface	177 x 63 cm
Height adjustment	46 cm - 93 cm
Backrest raises/ lowers	-20° to +85°
Footrest raises/ lowers	-85° to +15°
Electric max. load	260 kg
Trendelenburg	-15°

Electric Features:

Power supply	110-240v / 50-60
Maximum current intake	5A, 100-240V
Protection indicator	IPX4
Protection class	Class 2.B
Electric shock protection	TYPE B

Optional Accessories:

- « Extra Wide Upholstery (75 cm)
- « Breathing Hole
- « Foot Protection Cover
- « Detachable Head Rest
- « Electric Backrest
- « IV Pole attachment
- « Battery Operating System
- « Clean Edge Cover



						
BEIGE	CITRUS GREEN	APPLE MINT	CANDY	RASPBERRY CRUSH	DAMSON	PILLAR BOX RED
						
WHITE	ASH GREY	GUN METAL GREY	BLACK	DARK BLUE	MID BLUE	SKY BLUE

Medical Chairs

MED04EC, MED04ECT



Structure

Epoxy coated steel frame.
Electric backrest and gas foot rest as standard.
Industry-leading Linak 'High Reliability' actuators.
Electric lifting system provides optimum stability and supports weights of up to 190 kg.
Choice of tilting or non-tilting seat available.
Compact design for when space is limited.



Intended Use

The Medi-Plinth Single Column Medical Chair is a space-saving solution designed for a variety of clinical applications. Its innovative single-column design provides exceptional stability while allowing smooth height adjustments and effortless positioning.

This chair offers superior comfort for patients and ergonomic benefits for clinicians, making it ideal for procedures such as phlebotomy, examinations, and minor treatments. Durable upholstery, a robust frame, and optional accessories ensure it meets the demands of any medical environment.



Safety

MRSA resistant anti-bacterial & fire retardant upholstery.
Comprehensive parts and labour warranty covering manufacturing or material defects of seven years for the frame and three years for moving parts.
Complies with all applicable medical regulations.



Easy Cleaning

Upholstery compatible with all NHS approved cleaning solutions such as Chlor-Clean, Haz-Tab, Milton and Tristel, plus Clinell wipes and other alcohol based wipes.



Comfort

Soft, durable vinyl to ensure patient comfort. Its sturdy top makes the couch safe for the therapist to treat the patient with ease and ensures patient comfort.

Specifications

Single Column Medical Chairs

COLUMN MEDICAL CHAIRS

Column lift Medical Chair - Non-Tilting: **MED04EC**

Column lift Medical Chair - Tilting: **MED04ECT**

Patient surface	169 x 63 cm
Height adjustment	46 cm - 86 cm
Backrest raises	85°
Footrest lowers	- 85°
Backrest lowers (Tilting version)	-15°
Trendelenburg (Tilting version)	-15°
Electric max. load	190 kg

Electric Features:

Power supply	110-240v / 50-60
Maximum current intake	5A, 100-240V
Protection indicator	IPX4
Protection class	Class 2.B
Electric shock protection	TYPE B

Optional Accessories:

- « Extra Wide Upholstery (75 cm)
- « Foot Protection Cover
- « Foot Switch
- « Plaster Tables
- « Phlebotomy Arms
- « Adjustable Arm Rests



Medical Chairs

MED06E/H, MEDT06E



Structure

Epoxy coated steel frame.
Gas spring assisted back and leg sections.
Industry-leading Linak 'High Reliability' actuators.
Electric lifting system provides optimum stability and supports weights of up to 260 kg.
Hydraulic lifting system provides optimum stability and supports weights of up to 225 kg.
Choice of tilting or non-tilting seat available.



Intended Use

The Medi-Plinth General Medical Chair is a dependable and adaptable solution for a wide range of medical environments. Designed with both patient comfort and clinician efficiency in mind, this chair is perfect for procedures such as examinations, consultations, and minor treatments.

Featuring adjustable height, tilt, and backrest options, the chair ensures optimal patient positioning and ease of use for practitioners. Its sturdy construction, durable upholstery, and optional mobility features make it a practical choice for busy healthcare settings.



Safety

Anti-entrapment Suspa gas strut operated back and footrest.
MRSA resistant anti-bacterial & fire retardant upholstery.
Comprehensive parts and labour warranty covering manufacturing or material defects of seven years for the frame and three years for moving parts.
Complies with all applicable medical regulations.



Easy Cleaning

Upholstery compatible with all NHS approved cleaning solutions such as Chlor-Clean, Haz-Tab, Milton and Tristel, plus Clinell wipes and other alcohol based wipes.



Comfort

Soft, durable vinyl to ensure patient comfort. Its sturdy top makes the couch safe for the therapist to treat the patient with ease and ensures patient comfort.

Specifications

General Medical Chairs

GENERAL MEDICAL CHAIRS

Medical Chair Hydraulic - Non-Tilting: **MED06H**

Medical Chair Electric - Non-Tilting: **MED06E**

Medical Chair - Tilting: **MEDT06E**

Patient surface	186 x 63 cm
Height adjustment	46 cm - 91cm
Backrest raises/ lowers	-20 - 85°
Footrest lowers	-85°
Footrest raises (Tilting versions only)	15°
Hydraulic max. load	225 kg
Electric max. load	260 kg

Electric Features:

Power supply	110-240v / 50-60
Maximum current intake	5A, 100-240V
Protection indicator	IPX4
Protection class	Class 2.B
Electric shock protection	TYPE B

Optional Accessories:

- « Extra Wide Upholstery (75 cm)
- « Breathing Hole
- « Foot Protection Cover
- « Detachable Head Rest
- « Electric Backrest
- « Tilting Seat
- « Foot Switch
- « Daycare Armrests
- « Clean Edge Cover



Rehabilitation

T01E, T02E, T03E



Structure

Epoxy coated steel frame.
Gas spring assisted back and leg sections.
Industry-leading Linak 'High Reliability' actuators.
Electric lifting system provides optimum stability and supports weights of up to 180 kg.
Inclinometer included as standard.
3 Velcro harnesses and a single-piece fixed foot board included as standard.
Available as split leg, standard or two section tables.



Intended Use

Medi-Plinth's Tilt Tables are expertly engineered to support a wide range of clinical applications, including rehabilitation, diagnostics, and therapy.

Key features include a robust frame, durable upholstery, and secure patient straps for added safety during procedures. The addition of transport castors ensures effortless mobility, making these tables adaptable to different clinical spaces.



Safety

Anti-entrapment Suspa gas strut operated backrest (two sections models only).
MRSA resistant anti-bacterial & fire retardant upholstery.
Comprehensive parts and labour warranty covering manufacturing or material defects of seven years for the frame and three years for moving parts.
Complies with all applicable medical regulations.



Easy Cleaning

Upholstery compatible with all NHS approved cleaning solutions such as Chlor-Clean, Haz-Tab, Milton and Tristel, plus Clinell wipes and other alcohol based wipes.



Comfort

Soft, durable vinyl to ensure patient comfort. Its sturdy top makes the couch safe for the therapist to treat the patient with ease and ensures patient comfort.

Specifications

Tilt Tables

TILT TABLES

Split Leg Tilt Table - Electric: **T02E**

Tilt Table - Electric: **T01E**

2 Section Tilt Table - Electric: **T03E**



Patient surface	180 x 63 cm
Height adjustment	61.5 cm - 106.5 cm
Backrest raises/ lowers (Two Section Only)	-15 - 85°
Tilt adjustment	90°
Electric max. load	180 kg

Electric Features:

Power supply	110-240v / 50-60
Maximum current intake	5A, 100-240V
Protection indicator	IPX4
Protection class	Class 2.B
Electric shock protection	TYPE B

Optional Accessories:

- « Footswitch
- « Extra Wide Upholstery (75 cm)
- « Breathing Hole
- « Detachable Head Rest
- « Electric Backrest (2 Section Only)
- « Tilting Foot Board
- « Tilting Split Foot Board
- « Tilting Patient Table



						
BEIGE	CITRUS GREEN	APPLE MINT	CANDY	RASPBERRY CRUSH	DAMSON	PILLAR BOX RED
						
WHITE	ASH GREY	GUN METAL GREY	BLACK	DARK BLUE	MID BLUE	SKY BLUE

Rehabilitation

PHYS+02E/H, PHYS+03E/H



Structure

Epoxy coated steel frame.
Gas spring assisted back and foot section
Industry-leading Linak 'High Reliability' actuators.
Electric lifting system provides optimum stability and supports weights of up to 260 kg.
Hydraulic lifting system provides optimum stability and supports weights of up to 225 kg.
Single foot switch as standard.
Paper roll holder as standard.
Stitched key hole shaped breathing hole and bung.



Intended Use

These versatile plinths provide flexible positioning options, ensuring optimal comfort and support during physiotherapy, manual therapy, and a wide range of treatment sessions.

Ideal for physiotherapists, rehabilitation centers, and wellness clinics, Medi-Plinth's 2 & 3 Section Physio Plinths combine practical design with exceptional patient care, making them the perfect choice for modern treatment environments.



Safety

Anti-entrapment Suspa gas strut operated backrest/footrest.
MRSA resistant anti-bacterial & fire retardant upholstery.
Comprehensive parts and labour warranty covering manufacturing or material defects of seven years for the frame and three years for moving parts.
Complies with all applicable medical regulations.



Easy Cleaning

Upholstery compatible with all NHS approved cleaning solutions such as Chlor-Clean, Haz-Tab, Milton and Tristel, plus Clinell wipes and other alcohol based wipes.



Comfort

Soft, durable vinyl to ensure patient comfort. Its sturdy top makes the couch safe for the therapist to treat the patient with ease and ensures patient comfort.

Specifications

2 & 3 Section Physio Plinths

TWO SECTION

2 Section Physio Plinth - Electric: **PHYS+02E**

2 Section Physio Plinth - Hydraulic: **PHYS+02H**

Patient surface	186 x 63 cm
Height adjustment	46 cm - 91 cm
Backrest raises/lowers	-20° to +80°
Electric max. load	260 kg
Hydraulic max. load	225 kg

Electric Features:

Power supply	110-240v / 50-60
Maximum current intake	5A, 100-240V
Protection indicator	IPX4
Protection class	Class 2.B
Electric shock protection	TYPE B



THREE SECTION

3 Section Physio Plinth - Electric: **PHYS+03E**

3 Section Physio Plinth - Hydraulic: **PHYS+03H**

Patient surface	186 x 63 cm
Height adjustment	46 cm - 91 cm
Backrest raises/lowers	-20° to +80°
Electric max. load	260 kg
Hydraulic max. load	225 kg
Footrest raises/lowers	80°

Optional Accessories:

- « Handset
- « Extra Wide Upholstery (75 cm)
- « Foot Protection Cover
- « Cotsides
- « Adjustable Arm Rests



Rehabilitation

PHYS+04E/H, PHYS+05E/H



Structure

Epoxy coated steel frame.
Gas spring assisted back and foot section.
Electric lifting system provides optimum stability and supports weights of up to 260 kg.
Hydraulic lifting system provides optimum stability and supports weights of up to 225 kg.
Long stitched breathing hole.
Upholstery gusset linking all upholstery sections and covering all hinge points.
Single foot switch as standard.
Paper roll holder as standard.



Intended Use

Medi-Plinth's 4 and 5 Section Physio Plinths are designed to meet the demands of advanced physiotherapy and rehabilitation treatments. With multiple adjustable sections, these plinths offer unparalleled flexibility, allowing precise patient positioning for a variety of therapeutic and diagnostic procedures.

The additional sections provide enhanced support for complex positioning needs, making them ideal for specialised therapies.



Safety

Anti-entrapment Suspa gas strut operated backrest/footrest.
MRSA resistant anti-bacterial & fire retardant upholstery.
Comprehensive parts and labour warranty covering manufacturing or material defects of seven years for the frame and three years for moving parts.
Complies with all applicable medical regulations.



Easy Cleaning

Upholstery compatible with all NHS approved cleaning solutions such as Chlor-Clean, Haz-Tab, Milton and Tristel, plus Clinell wipes and other alcohol based wipes.



Comfort

Soft, durable vinyl to ensure patient comfort. Its sturdy top makes the couch safe for the therapist to treat the patient with ease and ensures patient comfort.

Specifications

4 & 5 Section Physio Plinths

FOUR SECTION

4 Section Physio Plinth - Electric: **PHYS+04E**

4 Section Physio Plinth - Hydraulic: **PHYS+04H**

Patient surface	186 x 63 cm
Height adjustment	46 cm - 91 cm
Backrest raises/lowers	-20° to +80°
Electric max. load	260 kg
Hydraulic max. load	225 kg

Electric Features:

Power supply	110-240v / 50-60
Maximum current intake	5A, 100-240V
Protection indicator	IPX4
Protection class	Class 2.B

FIVE SECTION

5 Section Physio Plinth - Electric: **PHYS+05E**

5 Section Physio Plinth - Hydraulic: **PHYS+05H**

Patient surface	186 x 63 cm
Height adjustment	46 cm - 91cm
Backrest raises/lowers	-20° to +80°
Electric max. load	260 kg
Hydraulic max. load	225 kg
Footrest raises	80°

Optional Accessories:

- « Handset
- « Extra Wide Upholstery (75 cm)
- « Foot Protection Cover
- « Cotsides
- « Positioning Wedge Set



Rehabilitation

PHYS+BOE



Structure

Epoxy coated steel frame.
Electric lifting system provides optimum stability and supports weights of up to 320 kg.
Stitched breathing hole and bung for patient.
Upholstery gusset linking all upholstery sections and covering all hinge points.
Single foot switch as standard.
Paper roll holder as standard.



Intended Use

Designed with simplicity and functionality in mind, this plinth provides essential features to support both patients and clinicians during therapy sessions.

It utilises the highest quality materials and innovative design to meet the diverse needs of physiotherapy practices. Offering great clinician access and positioning for larger patients, it provides a reliable platform for various physiotherapy exercises and treatments. NHS approved and infection control friendly, elevate your practice optimizing treatment outcomes and patient satisfaction.



Safety

MRSA resistant anti-bacterial & fire retardant upholstery.
Extra low entry height for wheel chair transfer.
Comprehensive parts and labour warranty covering manufacturing or material defects of seven years for the frame and three years for moving parts.
Complies with all applicable medical regulations.



Easy Cleaning

Upholstery compatible with all NHS approved cleaning solutions such as Chlor-Clean, Haz-Tab, Milton and Tristel, plus Clinell wipes and other alcohol based wipes.



Comfort

Soft, durable vinyl to ensure patient comfort. Its sturdy, wide top makes the couch safe for the therapist to treat the patient with ease and ensures patient comfort.

Rehabilitation

2 Section Physio Bo Plinth

PHYSIO BO

2 Section Physio Bo Plinth - Electric: **PHYS+BOE**

Patient surface	180 x 102 cm
Height adjustment	46 cm - 91 cm
Backrest raises/ lowers (Two Section Only)	-20° to +80°
Electric max. load	320 kg

Electric Features:

Power supply	110-240v / 50-60
Maximum current intake	5A, 100-240V
Protection indicator	IPX4
Protection class	Class 2.B
Electric shock protection	TYPE B

Optional Accessories:

- « Handset
- « Extra Wide Upholstery (120 cm)
- « Foot Protection Cover
- « Detachable Head Rest
- « Cotsides
- « Positioning Wedge Set



Rehabilitation

PHYS+MTE



Structure

Epoxy coated steel frame.
Electric lifting system provides optimum stability and supports weights of up to 320 kg.
Heavy duty adjustable foot for uneven floors.
Single foot switch as standard.
Paper roll holder as standard.



Intended Use

Medi-Plinth's Physio Mat Table is purposefully designed to provide a stable and spacious surface for physiotherapy and rehabilitation exercises. Its robust construction and generous dimensions make it ideal for treatments requiring full-body support, such as stretching, manual therapy, and mobility exercises.

With a strong steel frame and durable, easy-to-clean upholstery, this mat table ensures maximum patient comfort and long-lasting performance in busy clinical environments. The low height allows easy access for patients, while the table's sturdy design guarantees reliable stability for both therapist and patient during use.



Safety

MRSA resistant anti-bacterial & fire retardant upholstery.
Extra low entry height for wheel chair transfer.
Comprehensive parts and labour warranty covering manufacturing or material defects of seven years for the frame and three years for moving parts.
Complies with all applicable medical regulations.



Easy Cleaning

Upholstery compatible with all NHS approved cleaning solutions such as Chlor-Clean, Haz-Tab, Milton and Tristel, plus Clinell wipes and other alcohol based wipes.



Comfort

Soft, durable vinyl to ensure patient comfort. Its sturdy, wide top makes the couch safe for the therapist to treat the patient with ease and ensures patient comfort.

Specifications

Physio Mat Table

PHYSIO MAT TABLE

Mat Table Physio Bo Plinth - Electric : **PHYS+MTE**

Electric Features:

Patient surface	180 x 102 cm
Height adjustment	46 cm - 91 cm
Electric max. load	320 kg

Power supply	110-240v / 50-60
Maximum current intake	5A, 100-240V
Protection indicator	IPX4
Protection class	Class 2.B
Electric shock protection	TYPE B

Optional Accessories:

- « Handset
- « Extra Wide Upholstery (120 cm)
- « Cotsides
- « Positioning Wedge Set



Imaging

US03TE



Structure

Epoxy coated steel frame.
Industry-leading Linak 'High Reliability' actuators.
Stainless steel tray included as standard.
Electric lifting system provides optimum stability and supports weights of up to 260 kg.
The electric seat enables Trendelenburg.
Optional electric backrest available.



Easy Cleaning

Upholstery compatible with all NHS approved cleaning solutions such as Chlor-Clean, Haz-Tab, Milton and Tristel, plus Clinell wipes and other alcohol based wipes.



Comfort

Soft, durable vinyl to ensure patient comfort.
Patient-centric design enables enhanced positioning for improved comfort during examinations.



Safety

Anti-entrapment Suspa gas strut operated backrest/footrest.
MRSA resistant anti-bacterial & fire retardant upholstery.
Comprehensive parts and labour warranty covering manufacturing or material defects of seven years for the frame and three years for moving parts.
Complies with all applicable medical regulations.



Intended Use

The Medi-Plinth Ultrasound Couch is specifically designed to provide optimal support and comfort during diagnostic imaging procedures. With a sturdy frame, smooth height adjustment, and durable upholstery, this couch ensures patient stability and ease of access for clinicians during ultrasound examinations.

Ideal for use in hospitals, clinics, and diagnostic centers, the Ultrasound Couch combines reliability and functionality for high-quality patient care.

Specifications

ULTRASOUND COUCH

ULTRASOUND PLINTH

Ultrasound Plinth Electric - Tilting: **US03TE**

Patient surface	177 x 63 cm
Height adjustment	46 cm - 91 cm
Backrest raises	85°
Legrest lowers	-85°
Trendelenburg	-15°
Electric max load	260kg

Electric Features:

Power supply	110-240v / 50-60 Hz
Maximum current intake	5A, 100-240V
Protection indicator	IPX4
Protection class	Class 2.B
Electric shock protection	TYPE B

Optional Accessories:

- « Footswitch
- « Extra Wide Upholstery (75 cm)
- « Breathing Hole
- « Foot Protection Cover
- « Detachable Head Rest
- « Electric Backrest
- « Paper Roll Holder
- « Adjustable Arm Rests
- « Clean Edge Cover



						
BEIGE	CITRUS GREEN	APPLE MINT	CANDY	RASPBERRY CRUSH	DAMSON	PILLAR BOX RED
						
WHITE	ASH GREY	GUN METAL GREY	BLACK	DARK BLUE	MID BLUE	SKY BLUE

Imaging

CTABLE001



Structure

Fully-locking medical grade anti-static castors.
Carbon Fibre Top 208cm x 61cm (Length x Width).
Narrow base frame for ease of patient transfer.
Easy access emergency stop button for additional patient safety.
150cm carbon fibre imaging area.
Can be fitted with a wipeable, sealed anti-static mattress held in place by Velcro strips.
Manufactured in the UK and CE Marked.

The C Arm Imaging Table is available with a specific range of accessories that can provide positioning solutions for most procedures.



Easy Cleaning

Upholstery compatible with all NHS approved cleaning solutions such as Chlor-Clean, Haz-Tab, Milton and Tristel, plus Clinell wipes and other alcohol based wipes.



Safety

Linked brake castor system and kick brake castors.
Additional safety protection during procedures; the Stop Button can be used to isolate the lift mechanism.
Comprehensive parts and labour warranty covering manufacturing or material defects of seven years for the frame and three years for moving parts.



Intended Use

The Medi-Plinth C-Arm Imaging Table is designed to provide the ultimate support and stability during imaging procedures. With a robust, low-profile frame and exceptional adjustability, this table ensures optimal patient positioning and uninterrupted access for accurate imaging in clinical and surgical settings.

Its reliable design accommodates both general and specialized imaging techniques, making it an essential tool for hospitals, diagnostic centers, and surgical suites.

Specifications

C-ARM IMAGING TABLE

C-ARM IMAGING TABLE

C-Table Imaging Table: **CTABLE001**

Mattress density	75 cm
Patient surface	208 x 61 cm
Height adjustment	50 cm - 95 cm
Artefact free imaging area	150 cm
Electric max load	200 kg

Electric Features:

Power supply	110-240v / 50-60 Hz
Maximum current intake	5A, 100-240V
Protection indicator	IPX4
Protection class	Class 2.B
Electric shock protection	TYPE B
Backup lithium battery output	24V DC – 2.2 Ah

Optional Accessories:

- « 2x Full Length Accessory Rails
- « Lithium Battery Operating System
- « Stainless Steel Handle Bar
- « Mattress with Face Hole (Visco Foam/Anti-Static)



Please note mattress is not included as standard.



Fixed Height Plinths

PP02



Structure

Supports weights of up to 180kg.
Paper roll holder fitted as standard.
Epoxy coated steel frame.
Fixed height plinth, supplied flat pack.
Gusset between upholstery sections included to maintain effective infection control.
Available in a choice of 14 vinyl colours



Easy Cleaning

Upholstery compatible with all NHS approved cleaning solutions such as Chlor-Clean, Haz-Tab, Milton and Tristel, plus Clinell wipes and other alcohol based wipes.



Comfort

Soft, durable vinyl to ensure patient comfort. Patient-centric design enables enhanced positioning for improved comfort during examinations.



Safety

MRSA resistant anti-bacterial & fire retardant upholstery.
Comprehensive parts and labour warranty covering manufacturing or material defects of seven years for the frame and three years for moving parts.
Complies with all applicable medical regulations.



Intended Use

This plinth is a robust and cost-effective solution designed for everyday use in general practice, outpatient clinics, and treatment rooms. With its sturdy fixed-height frame and durable construction, this couch provides a stable and reliable platform for a variety of examinations and procedures.

Key features include high-quality, easy-to-clean upholstery for enhanced hygiene and patient comfort, as well as a minimalist design that prioritizes practicality and efficiency.

Specifications

PRACTICE PLINTHS


PRACTICE PLINTH

Two Section Practice Plinth: **PP02**


Patient surface	190 x 63 cm
Backrest raises	55°
Safe working load	180 kg



Top Elevation:



Operational Parameters:



UPHOLSTERY

All Medi-Plinth couches are available in our 14 standard colours. All products are upholstered in durable and comfortable vinyl which inhibits MRSA (3352 Methicillin resistant to Staphylococcus Aurous) Bio-Pruf treated to resist the growth of mould, mildew and odour causing microbes.

Guide colours are for indication use only. Samples available on request.



Apple Mint



Citrus Green




Sky Blue



Mid Blue




Dark Blue



Damson




Raspberry Crush



Pillar Box Red



Candy




Beige



Ash Grey



Gun Metal Grey



Black



White



Medi Plinth

www.medi-plinth.co.uk

EU-MDR | UK
CA | CE

7-11 Holywells Road, Ipswich, Suffolk IP3 0DL - Call: +44 (0) 01473 212 010