

Flexline Phlebotomy Cart



FLXIV Phlebotomy Cart

The Metro Flexline Phlebotomy Cart is built to enhance efficiency, security, and hygiene within healthcare-grade blood collection workflows.

Built from lightweight polymer with ergonomic push handles and smooth-rolling casters, it moves easily through busy healthcare environments and positions neatly at the point of care. The self-closing drawers offer full access and can be secured with serialized plastic security tags to maintain chain-of-custody integrity, while both a key lock and an electronic keypad provide controlled access for authorised personnel.

The cart features a practical combination of three drawer depths (7.5 cm, 15.2 cm, and 22.9 cm) paired with dividers, an ampule insert tray, an IV pole, and a sharps bin holder to keep essential items organised and within reach. A generous flat work surface on top supports preparation and procedural tasks, while Microban antimicrobial protection helps maintain infection control. The cart balances compact form with sufficient capacity, making it an efficient, durable, and secure solution for phlebotomy and other point-of-care applications.

Code	Width	Length	Height	Weight
FLXIV	82cm	56cm	98cm	80kg

Configuration

FLXIV: Phlebotomy Cart with Key Lock System & Electronic Keypad Locking System, 3 x 7.5Cm, 1 x 15.2cm & 1 x 22.9cm Deep Drawers with Red Drawer Pulls. Accessories: 2 x 75mm Drawer Divider Kit, 1 x 152/229mm Drawer Divider Kit, 1 x Ampule Insert & Divider Tray, IV Pole & Sharps Bin Holder

Accessory	Code
3" Drawer Divider Kit for Flexline	FL113
6"/9" Drawer Divider Kit for Flexline and Lifeline Carts	FL116
Flexline Oxygen Tank Holder, All Size Carts, 4" Diameter Tank	FL234
Adjustable Defibrillator Tray for Lifeline and Flexline Carts	LEC304