

# Medical Recliner Chair



## 6801 Hydraulic Recliner Chair

Crafted for medical examination, sampling, and safe patient transfers, the Hydraulic Medical Recliner Chair incorporates a stepless hydraulic height adjustment from 51cm to 83cm, operated via a convenient foot pump.

The backrest and leg section are individually adjustable into a fully reclined position with gas spring assistance, while a folding footrest ensures both patient comfort and space efficiency.

Its robust epoxy-coated steel frame and seamless, hygienic upholstery meet stringent infection-control requirements, and the chair includes features like locking 125mm castors with directional control, five-position folding armrests with quick-lock means, an ergonomically shaped seat, and a neck cushion with hygiene-oriented handles and knobs.

With compact maximum dimensions—67cm width and 159–168cm length this chair excels in both small spaces and demanding clinical environments. Its mobility is enhanced by durable castors and optional front wheels, while accessory integration, including IV pole adapters and a hygienic sampling pillow, ensure comprehensive functionality.

Code	Total Width	Total Length	Seat Width	Seat Depth	Backrest Adj
6801 AM	67cm	159 - 168cm	47cm	46cm	90°

Code	Leg Section Adj	Headrest Height Adj	SWL	Weight w/o Accessories
6801 AM	90°	9cm	150kg	72kg

### Configuration

6801 AM: Medical Recliner Chair, Hydraulic- inc Neck Cushion, Meal Tray, Adjustable Length Foot Rest, Arm Rests & Push Handles

Accessories	Code
Ø300mm Wheels in Front	WLFW300/125GH
Calf Support Length Adjustment	Calf support (slide)
IV Pole Adapter	GH20871
IV Pole, 4 Hooks. Height Adjustment with One Hand - Adapter Required	60121AMI
Durable Hygiene Pillow for Sampling 35x45cm	50045