

Starsys Catheter Storage Cabinet



SXRS76MCATH1/ SXR76MCATH1

The Starsys Mobile Catheter Storage Cabinet is designed to deliver an efficient, hygienic, and highly adaptable solution for catheter management in healthcare environments. Its streamlined, modular system combines aesthetics, cleanliness, and operational flexibility which is ideal for busy wards, theatres, or transport hubs.

Constructed from advanced polymer with integrated Microban® antimicrobial protection, these units resist bacterial growth while simplifying routine cleaning and disinfection protocols. Standard configurations include full-extension catheter slides with visual labels and sturdy hooks, enabling quick access and visual inventory control. The single-unit model houses one shelf, three slides, 15 hooks, and 30 labels, supporting storage for up to 180 catheters, with capacity expandable using additional hooks. The double-unit version doubles the system, offering two shelves, six slides, 30 hooks, and space for 360 catheters.

Cabinets include secure left-hinged 270° swing doors, each with an M300-series key lock and two keys, ensuring safe storage during movement and transit.



Code	Width	Depth	Height	Weight
SXRS76MCATH1	57.6cm	63.2cm	199cm	100kg
SXR76MCATH1	107.9cm	63.2cm	199cm	136kg

Configuration

SXRS76MCATH1: Mobile Single Catheter Unit with 1 Catheter Shelf, Locking Doors, 3 Slides, 15 Hooks and 30 Labels. Capacity 180 Catheters

SXR76MCATH1: Mobile Double Catheter Unit with 2 Catheter Shelf, Locking Doors, 6 Slides, 30 Hooks and 30 Labels. Capacity 360 Catheters

