

Lifeline® Emergency Carts



LECCRP4 Emergency Carts

The Metro LECCRP4 Lifeline Crash Cart is purpose-built for fast, reliable response in emergency care environments. Its advanced polymer construction with integrated Microban® antimicrobial protection ensures long-lasting durability and supports strict infection control standards. Integrated push handles, and large 5-inch swivel casters with a fifth-wheel steering assist provide effortless manoeuvrability in tight or busy clinical spaces.

The cart features two 3-inch, one 6-inch, and one 12-inch deep drawer, allowing rapid access to essential equipment. Tamper-evident security seals help prevent unauthorised access.

Each unit comes complete with side storage bins, a defibrillator tray, suction pump shelf, IV pole mount, oxygen tank holder, and backboard kit.

Available in both blue and red drawer pull options, the LECCRP4 allows departments to colour-code carts for specific clinical functions, improving workflow and reducing critical response time.



Code	Width	Length	Height	Weight
LECCRP4	57.4cm	96.5cm	106.9cm	90.7kg

Configuration

LECCRP4XP: Lifeline Crash Cart with Passive Lock System, 2 x 7.5cm, 1 x 15.2cm & 1 x 30.5cm Deep Drawers with Blue Drawer Pulls. Inc. 3 Locking Side Bins, Defib Platform, Backboard, IV Pole, Gel Bin, Security Seals

LECCRP4XPR: Lifeline Crash Cart with Passive Lock System, 2 x 7.5cm, 1 x 15.2cm & 1 x 30.5cm Deep Drawers with Red Drawer Pulls. Inc. 3 Locking Side Bins, Defib Platform, Backboard, IV Pole, Gel, Security Seals

Accessory	Code
3" Drawer Divider Kit for Flexline Lifeline Carts	FL113
6"/9" Drawer Divider Kit for Flexline and Lifeline Carts	FL116
Top Cavity Tray	LEC143
Plastic Security Seals for Lifeline Carts, 100 Pack	LEC320
Suction Pump Shelf for Lifeline and Flexline Carts	LEC306