

Rollators



Polo Plus & Polo Comfort Rollators

The Polo Plus-T and Polo-Comfort rollators are high-quality mobility solutions designed to combine safety, durability, and ease of use in both healthcare and community settings.

The Polo Plus-T Rollator is a lightweight aluminium folding rollator engineered for safety, durability, and ease of use in both healthcare and community environments. It incorporates a dual-action braking system, puncture-proof tyres, and a reversible resting seat, making it ideal for short breaks during use. An integrated Velcro-secured storage pouch offers convenient space for personal items. Available in small, regular, and tall sizes, the Polo Plus-T accommodates a range of user heights while maintaining a stable footprint.

The Polo-Comfort Rollator builds on the robust design of the Polo Plus-T while enhancing user experience with a padded seat and backrest, offering additional comfort for prolonged use. Retaining the same folding frame, braking reliability, and ergonomic sizing, it provides the same capacity and stability while offering improved seated support.



Code	Width	Depth	Height of Handle	Weight	SWL	Seat Height	Seat Width
181.10	67cm	70cm	64-74cm	7.4kg	120kg	49cm	48cm
181.20	67cm	70cm	69-82cm	7.6kg	120kg	55cm	48cm
181.30	67cm	70cm	74-92cm	7.8kg	120kg	60cm	48cm
182.20	68cm	73cm	70-83cm	7.8kg	120kg	56cm	48cm
182.30	68cm	78cm	77-96cm	8kg	120kg	62cm	48cm

Configuration

- 181.10: Polo Plus-T Aluminium Rollator with Bag, Size S
- 181.20: Polo Plus-T Aluminium Rollator with Bag, Size M
- 181.30: Polo Plus-T Aluminium Rollator with Bag, Size L
- 182.20: Polo Comfort Aluminium Rollator Fully Equipped, Size M
- 182.30: Polo Comfort Aluminium Rollator Fully Equipped, Size L



Find our products on the NHS Supply Chain Catalogue

