

Ward Case Note Trolleys



CNT008/09/10/11

The Shuttleworth Case Notes Trolley offers a secure, hygienic, and highly manoeuvrable solution for managing patient notes and medical imaging within healthcare environments. Built on a robust, powder-coated mild steel frame with aluminium top, the trolley features solid side panels and a lower storage shelf, combining strength with a streamlined design for ease of cleaning and infection control.

The top section of the trolley includes a lockable cabinet, available with either Digital or Cam lock options to safeguard sensitive documents and comply with patient confidentiality standards.

The trolley's mobility is enhanced by four 100mm swivel castors, complete with two braked for stability, and ergonomically placed side push handles ensure smooth navigation and precise positioning, even in crowded corridors or tight spaces.

Making it ideal for ward, clinic, or administrative settings where efficiency and security are essential.

Code	Width	Depth	Height
CNT008	80cm	46.7cm	94.8cm
CNT009	80cm	46.7cm	94.8cm
CNT011	100cm	46.7cm	94.8cm

Configuration

CNT008: Case Notes Trolley with Cam Lock

CNT009: Case Notes Trolley with Digital Lock

CNT011: Case Notes Trolley, Large with Digital Lock



Find our products on the NHS Supply Chain Catalogue

