

Emergo Patient Trolley



6280 & 6380 Patient Trolleys

The Emergo patient trolleys are purpose-built to deliver reliable performance in patient transfer, examination and treatment across multiple hospital departments.

Constructed with a strong steel frame and high-quality components, these models combine stability, safety and ease of handling to withstand the heavy daily use of acute healthcare environments. Each trolley features a smooth, ergonomic design for infection control, with drop-down side rails for patient safety and height-adjustable sections to accommodate clinical procedures.

6280 features a two-section X-ray-compatible platform with subframe, streamlining imaging workflows by allowing scans to be conducted without having to move the patient. The 6380 enhances versatility with a three-section X-ray platform plus subframe, enabling independent backrest and leg section adjustment, ideal for more complex clinical positioning needs.

Both models are equipped with centrally locking, 150 mm antistatic castors. Designed for long-term reliability, these trolleys support a range of accessories including IV poles, oxygen cylinder holders, baskets and trays, making them an adaptable investment for healthcare facilities.

Code	Width	Patient Surface Width	Length	Weight	Height Adj	SWL
6280	75cm	70.5cm	215cm	105kg	62.5-99.5cm	225kg
6380	75cm	70.5cm	215cm	110kg	62.5-99.5cm	225kg

Code	Trendelenburg	Anti-Trendelenburg	Back Adj	Leg Adj	Base Sections
6280	14°	7°	0 - 70°	N/a	2
6380	14°	7°	0 - 70°	0 - 45°	3



Configuration

6280: 2-part X-ray bed w/ subframe, hydraulic height adjustment, gas assisted Trendelenburg & back section

6380: 3-part X-ray bed with hydraulic height adjustment, gas assisted Trendelenburg, back section & leg section



Antistatic Trolley Mattress (Accessory)

Trolley mattresses are designed for transport and X-ray stretchers. The interior is made of fire-resistant flexible foam. The interior of model 51701 (10 cm) consists of two layers: on the surface there is a softly forming visco-elastic layer to equalize tissue pressure. The top layer is a PU-coated antistatic black hygienic fabric that is impermeable to liquids but breathable. The zipper is protected by a drip flap. The mattress has Velcro strips on the bottom and four Velcro strips on the sides to keep the product in place.

Accessories	Code
Wire Basket Under Back Section	10000582
Paper Roll Holder	10000733
Long Push Bars, Detachable (pair)	10006770
Patient Restraining Straps (set)	10000790
Push bar, Foldable	10006793
Oxygen Bottle Holder (5kg Bottles)	10000589
Oxygen Bottle Holder (2kg Bottles)	100005860
Infusion Rod, 4-hooks	10000499
Hygiene Antistatic Mattress (50mm, 80mm or 100mm thick)	100051751
Trolley Bed End	10000567
Wire Basket (requires trolley end)	128009784

